

# PRO1 MICRO

Rev. 15 / 03.2024


**Chemdye®**

## Hygiene Monitoring System



### Instructions for use


Shake downwards vigorously until the solution turns to green and reaches the readout cone reference volume

Slide swab upwards without completely removing

### Results reference guide

**NOTE:** Reference colors identified in printed prospectus and boxes as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

<b>CE</b>	CE mark.	<b>LOT</b>	Batch number.
<b>EC REP</b>	Authorized representative in the European Community.	<b>MFG</b>	Manufacture Date.
<b>MD</b>	Medical Device.	<b>EXP</b>	Expiration Date.
		<b>MFG</b>	Manufacturer.

## EN Hygiene Monitoring System

*For the detection and quantification of surface proteins*

#### Product description

Reusable medical devices are required to be safely reprocessed prior to patient use. This includes cleaning, disinfection and sterilization. Protein soils are the biggest contamination source on instruments surfaces and the cleaning step is crucial for the appropriate protein removal.

Residual protein removal is pertinent to reduce the risk of prion disease transmission. Prion proteins are resistant to several sterilization processes and are the transmission agent of iatrogenic spongiform encephalopathies such as the variant Creutzfeldt-Jakob disease (vCJD)<sup>1,2</sup>. Moreover, instrument’s soil can act as a protectant and allow the survival of virus and bacteria to disinfection and sterilization processes<sup>3</sup>. Soiled proteins are the greatest challenge to the cleaning procedure due to its high adhesion, especially when they undergo a denaturation process, for example in a thermo-disinfection process.

In virtue of these facts, a fast and sensitive protein detection method is indispensable to evaluate the initial cleaning process. Total residual protein level is an indicator of the presence of potential toxic iatrogenic proteins and, additionally, it can be used

to predict the efficiency of further disinfection and sterilization processes.

Chemdye® PRO1 MICRO Hygiene Monitoring System is the most suitable tool for the verification of medical device’s hygiene. Chemdye® PRO1 MICRO has been designed to detect and quantify residual proteins and iatrogenic agents on medical instrument’s surfaces after the cleaning process. The system includes a high absorption swab, which allows the collection of samples from different types of surfaces with high efficiency. A visual readout of a color change indicates the presence of detectable levels of protein. Chemdye® PRO1 MICRO has a quantification range from 1 to 50 µg of protein with a limit of detection (LOD) of 0.5 µg\* in only 4 minutes.

#### Characteristics

Chemdye® PRO1 MICRO along with Bionova® Auto-readers provide all the necessary equipment for the Residual Protein Test described in HTM 01-01 and HTM 01-05 UK guidelines, APSIC guideline (Asia Pacific Society of Infection Control) and ISO 15883-1, 15883-5 Standards.

#### Indications for use

Use for testing any surface after the cleaning process and for hard-to-reach areas in complex instruments, such as serrated edges or box joints.

#### Precautions

Do not open until use. Do not reuse. Do not touch the swab or the inside of the device with your hands. Wear gloves, surgical mask and other hygiene measures. Do not freeze. High concentrations of alkaline detergents may cause false negative results in some cases. Hydrogen peroxide-based disinfectants may cause a false positive color change reaction. Although Chemdye® PRO1 MICRO has high sensitivity, it is not a microbiological monitoring device, and the surfaces should not be considered sterile according to the results obtained.

#### Instructions for use

- Remove the device from the pouch and **allow it to reach a temperature of 20-25 °C**.
- Remove the swab from the device (Fig. A). Carefully apply 2 drops of moisturizer on the swab tip.
- Swab thoroughly the surface to be tested (Fig. B). Swab in zigzag in one direction, and then do the same procedure perpendicularly (Fig. B). Press the swab and rotate it several times, while collecting the sample. To perform comparable hygiene tests, standardize the surface for sampling. In order to achieve this, it is advisable to always take the sample from the same surface region of the device, especially in the most hard-to-reach areas.
- Place the swab back into the device (Fig. C).
- Activate it, by pressing it firmly downwards (Fig. D).
- Shake vigorously downwards, until the solution turns to green and reaches the reference volume of the readout cone (Fig. E).
- Shake vigorously downwards for 15 seconds, with the swab submerged in the solution (Fig. F).

**WARNING:** Once activated, incubate the pen immediately for final readout.
8. Slide the swab upwards without removing it from the device (Fig. G).

9. Shake vigorously downwards and collect the solution until reaching the volume level indicated in the reference readout cone (Fig. H and I).

**CAUTION:** If the device is not shaken properly (so as the solution in the internal compartment falls completely), or if there is a great presence of bubbles in the solution, the results can be invalid.

- Make sure that the swab is outside the readout cone (Fig. I).
- Incubate the Chemdye® PRO1 MICRO device at 60 ± 2 °C (Fig. J) for:

-4 minutes in Bionova® MiniPro incubator.

-7 minutes in Bionova® IC10/20, IC10/20FR, IC10/20FRLCD incubators.

**IMPORTANT:** The incubation must be done without the swab submerged in the solution and immediately after the activation stage.

12. For a quantitative analysis, you may use Bionova® IC10/20FR, IC10/20FRLCD or MiniPro Auto-reader incubators. Read the instructions for use before using these devices.

13. For a semi-quantitative analysis, remove the Chemdye® PRO1 MICRO immediately after the incubation time and interpret visually by using the reference color guide.

**IMPORTANT:** The estimations should be done before 5 minutes after removing the device out of the incubator. Results analyzed beyond this period are not valid.

14. After using the Chemdye® PRO1 MICRO System, reprocess the instrument or the surface used for protein monitoring.
**IMPORTANT:** During the activation of PRO1 MICRO Hygiene Monitoring System, the swab tip could become detached. If this happens, repeat the test using a new pen. If the swab tip is retained in the reading cone, a visual estimation can be carried out after incubation at 60 °C, as detailed in section 12.

**NOTE:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

#### Storage

Chemdye® PRO1 MICRO System should be stored at temperatures between 2-30 °C, away from light and relative humidity (RH) of 30-80 %.

Store vertically as indicated on the packaging. Shelf life of 24 months when stored refrigerated (between 2-8 °C). Maximum shelf life at room temperature is 6 months (up to 30 °C). 1 month of storage at room temperature reduces the expiration date by 4 months.

#### Disposal

Discard the product after use according to your country’s healthcare and safety regulations.

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#### References

**↑** Bonda DJ, Manjila S, Mehndiratta P, et al. Human prion diseases: surgical lessons learned from iatrogenic prion transmission.

*Neurosurg. Focus.* **2016**;41(1):E10.

**↑** Hervé RC, Hedges J, Keevil CW. Improved surveillance of surgical instruments reprocessing following the variant Creutzfeldt-Jakob disease crisis in England: findings from a three-year survey. *J. Hosp. Infect.* **2021**;110:15-25.

**↑** Alfa MJ, DeGagne P, Olson N, Hizon R. Comparison of liquid chemical sterilization with peracetic acid and ethylene oxide sterilization for long narrow lumens. *Am. J. Infect. Control.* **1998**;26(5):469-77.

<p>Learn more about this product and check the digital color reference guide</p>  <p>Scan the QR code</p>
<p><b>Terragene S.A</b> Ruta Nacional N° 9, Km 280 - CP 2130. Parque Industrial Micropi- Alvear-Santa Fe-Argentina.</p>
<p><b>Terragene®</b></p>
<p><b>MD</b> <b>CE</b></p>