Result reference guide

Unprocessed

WARNING SHEET


PRIMARY SHEET

Unsatisfactory

WARNING SHEET

PRIMARY SHEET

Satisfactory

WARNING SHEET

PRIMARY SHEET

⚠ Sterilizer may be used, but maintenance should be scheduled.

Explanation of Symbols

(SEAM) AIR REMOVAL

Product designed for testing air removal in dynamic-air-removal steam sterilizers.

CE CE mark.

EC REP Authorized representative in the European Community.

LOT Batch number.

MFG Manufacture Date.

EXP Expiration Date.

MFG Manufacturer.

EN Bowie-Dick Test Pack

For vacuum-assisted steam sterilizers

Indications for Use

Terragene® Chemdye® Bowie-Dick Test Pack BD125X/1 is designed for testing air removal efficiency for cycles of 132 °C for 4 minutes and 134 °C for 3.5 minutes in dynamic-air-removal steam sterilizers.

Product description

Chemdye® Bowie-Dick Test Pack is used to evaluate the effectiveness of air removal during pre-vacuum steam sterilization cycles operating at 134 °C, 3.5 minutes / 132 °C, 4 minutes.

The Chemdye® Bowie-Dick Test Pack consists of chemical indicators on two test sheets (125 x 113 mm) positioned between porous materials and a reticulated foam, wrapped all together forming a package.

The test pack has a Type 1 chemical indicator printed on label. The primary test sheet with a regular pattern corresponds to a Bowie-Dick air-removal indicator (Type 2 according to ISO 11140-1). It turns from yellow to uniform dark brown/black color when processed. An air removal failure is indicated as a lighter-colored area in the center of the sheet (a non-uniform color change) indicating sterilizer malfunction.

The Test Pack has also a warning sheet, which is placed at a lower point in the pack than the primary test sheet. The warning sheet has an indicator ink circle that turns from yellow to dark brown/black for early detection of air removal failures before they appear on the primary test sheet, therefore allowing rapid detection of failures in the sterilization cycle.

The Bowie-Dick Test Pack comply to ANSI/AAMI/ISO 11140-1:2014 and ANSI/AAMI/ISO 11140-5:2007 standards.

IMPORTANT: Bowie-Dick tests should be used daily, before running the first load of the day, after a sterilizer is installed or relocated, after a sterilizer malfunction, after sterilization process failures and after any major repairs of the sterilizer. A shortened cycle (i.e., a cycle omitting the post-vacuum drying phase) should be run first to properly heat up the sterilizer.

WARNING: extending the manufacturer's recommended time for the Bowie-Dick Test may result in false endpoints developing because entrapped air diffuses over time. Do not use a Bowie-Dick Test Pack in gravity displacement steam sterilization cycles (where there is no pre-vacuum phase to check) or in other kinds of sterilization processes other than vacuum-assisted.

Instructions for use

1. Run a warm-up cycle to properly heat the sterilizer.
2. Place the Bowie-Dick Test Pack horizontally with the label facing up in an empty vacuum-assisted sterilizer chamber, in the bottom/front of the sterilizer rack and over the drain. The Test Pack should read "THIS SIDE UP" when placed correctly.
3. Run the sterilizer cycle for 3.5 min. - 134 °C / 4 min. - 132 °C.
4. Remove the test pack from sterilizer and allow it to cool down. Check color change of type 1 indicator on the label.
5. Unwrapped the pack and remove the test and warning sheets. Check them for color change and record the results.

NOTE: if any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Interpretation of results

After completing the cycle, retrieve and examine the test sheets using the Result Reference Guide.

Satisfactory test: Test sheet shows uniform dark brown/black color change. This indicates adequate air removal from the sterilization chamber.

Unsatisfactory test: This failure is caused by air pockets, incomplete air removal from the chamber, poor quality steam penetration and/or non-condensable gases in steam supply. If test sheet shows non-uniform color change, it should be reported immediately and stop using the sterilizer until examination for possible malfunction. Recall and quarantine the previous loads sterilized in that machine until you are certain that the machine has not failed before.

Warning sheet

If the indicator on the warning sheet does not turn to dark brown/black as the reference color, potential sterilizer failure is indicated and routine sterilizer maintenance should be scheduled. Refer to Result reference guide.

Storage

Store in a dry place, keep away from sunlight, Temperature: 10-30 °C, relative humidity: 30-80 %.

Do not store near sterilizing agents. It is recommended to protect the test sheets after use inside a plastic foil for storage.

Shelf life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both, product expiration date and batch number are specified on product packaging and in each individual card. Do not use indicators after their expiration date.

Endpoint stability: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard indicators after use along with paper waste according to your country's healthcare regulations.

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ATTENTION: If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive Trazanto Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes.

This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through Bionova® Software. Use Trazanto scanner to interpret the results (see instructions for use in the manual of the device).

ATTENTION: Do not write on the indicator before scanning.

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NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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