
BD8948X/1


Rev. 1 / 07.2021

Bowie-Dick Test Card


Result reference guide

Unprocessed


Failures

Purple shades or heterogeneous color change to green due to air removal failure.

Processed and correct

PANTONE 346 C or similar

NOTE: Reference colors exhibited in printed prospects and boxes as well as those displayed on our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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Explanation of Symbols

STEAM AIR REMOVAL

Product designed for testing air removal in dynamic-air-removal steam sterilizers.

CE CE mark.

EC REP Authorized representative in the European Community.

Indicators 100 % Metal Free

LOT Batch number.

MANUFACTURE Date.

EXP Expiration Date.

MANUFACTURER Manufacturer.

EN

Bowie-Dick Test Card

Effectiveness of air removal in vacuum-assisted Steam sterilizers

Indications for Use

Terragene® Chemdyne® Bowie-Dick Test Card BD8948X/1 is designed for testing air removal efficiency for cycles of 132 °C for 4 minutes and 134 °C for 3.5 minutes in dynamic-air-removal steam sterilizers.

Device description

Chemdyne® Bowie-Dick Test Card is used to evaluate the effectiveness of air removal during pre-vacuum steam sterilization cycles operating at 134 °C, 3.5 minutes / 132 °C, 4 minutes.

Chemdyne® Bowie-Dick Test Kit contains a BD8948X/1 air-removal indicator card (Type 2 according to ISO 11140-1) to be used with its reusable holder. The indicator turns uniformly from purple to green when processed correctly. If a heterogeneous color change occurs in the indicator, it is an indication of the presence of residual air in the chamber and, therefore, sterilizer malfunction. Cycle information can be recorded on the back of the indicator, thus allowing a permanent record of the process.

The Bowie-Dick Test Card comply to ANSI/AAMI/ISO 11140-1:2014 and ANSI/AAMI/ISO 11140-5:2007 standards.

IMPORTANT: Bowie-Dick Tests should be used daily, before running the first load of the day, after a sterilizer is installed or relocated, after a sterilizer malfunction, after sterilization process failures and after any major repairs of the sterilizer. A shortened cycle (i.e., a cycle omitting the post-vacuum drying phase) should be run first to properly heat up the sterilizer.

WARNING: extending the manufacturer's recommended time for the Bowie-Dick Test may result in false endpoints developing because entrapped air diffuses over time. Do not use a Bowie-Dick Test Card in gravity displacement steam sterilization cycles (where there is no pre-vacuum phase to check) or in other kinds of sterilization processes other than vacuum-assisted.

Instructions for use

1. Run a warm-up cycle to properly heat the sterilizer.
2. Insert the Bowie-Dick Test Card in its holder or in a similar suitable one with the chemical indicator side facing up. Place it in horizontal position on the bottom rack, near the drain of the empty sterilizer.
3. Run the sterilizer cycle for 3.5 min. - 134 °C / 4 min. - 132 °C.
4. Remove the test card from sterilizer and check it for color change.
5. Record the results on the back of the card and keep it as a record.

NOTE: if any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Interpretation of results

After completing the cycle, retrieve and examine the test card using the Result Reference Guide.

Satisfactory test: Card shows uniform green color change. This indicates adequate air removal from the sterilization chamber.

Unsatisfactory test: Card shows heterogeneous color change (purple remains). This failure is caused by air pockets, incomplete air removal from the chamber, poor quality steam penetration and/or non-condensable gases in steam supply. If test card shows non-uniform color change, it should be reported immediately and stop using the sterilizer until examination for possible malfunction. Recall and quarantine the previous loads sterilized in that machine until you are certain that the machine has not failed before.

Storage

Store in a dry place, keep away from sunlight, Temperature: 10-30 °C, relative humidity: 30-80 %.

Do not store near sterilizing agents.

Shelf Life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both, product expiration date and batch number are specified on product packaging and in each individual card. Do not use indicators after their expiration date.

Endpoint Stability: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard indicators after use along with paper waste according to your country's healthcare regulations.

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