

Bowie-Dick Test Sheet

● Unprocessed

● Unsatisfactory

● Satisfactory

Indicators 100% Metal Free.

Indicators for use in Bowie and Dick-type test.

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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Explanation of Symbols

Product designed for testing air removal in dynamic-air-removal steam sterilizers.

CE mark.

Authorized representative in the European Community.

Batch number.

Manufacture Date.

Expiration Date.

Manufacturer.

Indications for Use

Terragene® Chemdye® Bowie-Dick Test Sheet BDA4/1 is designed for testing air removal efficiency for cycles of 132 °C for 4 minutes and 134 °C for 3.5 minutes in dynamic-air-removal steam sterilizers.

Device description

Chemdye® Bowie-Dick Test Sheet is sized A4 air-removal indicator (Type 2 according to ISO 11140-1) in a concentric pattern to be used to evaluate the effectiveness of air removal during pre-vacuum steam sterilization cycles operating at 134 °C, 3.5 minutes / 132 °C, 4 minutes.

The sheets are designed to be used in standard ANSI/AAMI ST79 Bowie-Dick test packs. Homogenous color change from yellow to uniform dark brown/black color when processed indicates adequate air-removal and steam penetration. An air removal failure is indicated as a lighter-colored area in the center of the sheet (a non-uniform color change) indicating sterilizer malfunction.

The Bowie-Dick Test Sheet comply to ANSI/AAMI/ISO 11140-1:2014 and ANSI/AAMI/ISO 11140-5:2007 standards.

IMPORTANT: Bowie-Dick Tests should be used daily, before running the first load of the day, after a sterilizer is installed or relocated, after a sterilizer malfunction, after sterilization process failures, and after any major repairs of the sterilizer. A shortened cycle (i.e., a cycle omitting the post-vacuum drying phase) should be run first to properly heat up the sterilizer.

WARNING: extending the manufacturer's recommended time for the Bowie-Dick Test may result in false endpoints developing because entrapped air diffuses over time. Do not use a Bowie-Dick Test Sheet in gravity displacement steam sterilization cycles (where there is no pre-vacuum phase to check) or in other kinds of sterilization processes other than vacuum-assisted.

Instructions for use

1. Run a warm-up cycle to properly heat the sterilizer
2. Prepare the test pack of cotton surgical towels freshly laundered, but not ironed. The package should be conditioned in accordance with technical guidelines (ANSI/AAMI ST79 or ISO 11140-5) with standardized dimensions: (250 ± 20) mm x (300 ± 20) mm.
3. Place a single Chemdye® BDA4/1 Bowie-Dick Test Sheet in the center of the test pack.
4. Secure the pack with tape not exceeding 25 mm in width.
5. Place the test pack horizontally in the lower shelf of an empty vacuum assisted steam sterilizer.
6. Run the sterilizer cycle for 3.5 min. - 134 °C / 4 min. - 132 °C.
7. After completion of the cycle, remove the test sheet from the pack and examine for compliance.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Interpretation of results

After completing the cycle, retrieve and examine the test sheets using the Result Reference Guide.

Satisfactory test: test sheet shows uniform dark brown/black color change. This indicates adequate air removal from the sterilization chamber.

Unsatisfactory test: this failure is caused by air pockets, incomplete air removal from the chamber, poor quality steam penetration and/or non-condensable gases in steam supply. If test sheet shows non-uniform color change, it should be reported immediately and stop using the sterilizer until examination for possible malfunction. Recall and quarantine the previous loads sterilized in that machine until you are certain that the machine has not failed before.

Storage

Store in a dry place, keep away from sunlight, temperature: 10-30 °C, relative humidity: 30-80 %.
Do not store near sterilizing agents.

Shelf Life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both, product expiration date and batch number are specified on product packaging and in each individual test sheet. Do not use indicators after their expiration date.

Endpoint Stability: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard indicators after use along with paper waste according to your country's healthcare regulations.

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