


## Bowie-Dick Test Sheet


## ● Unprocessed

## ● Unsatisfactory

## ● Satisfactory

 Indicators 100 % Lead Free

Indicators for use in Bowie-Dick type test.

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NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

## Explanation of Symbols

 Product designed for for use with Steam sterilization cycles.

 CE mark.

 Authorized representative in the European Community.

 Medical Device.

 Batch number.

 Manufacture Date.

 Expiration Date.

 Manufacturer.

 **Bowie-Dick Test Sheet**

## Product description

Chemdye® BDA4 Bowie-Dick Test Sheets are A4 sized single-sheet printed with a Type 2 Chemical Indicator, in a concentric pattern. This type of sheet is intended to monitor the sterilizer, where the

main goal is to detect incorrect air removal, inadequate steam penetration and/or vacuum bomb failures of the said equipment.

A homogeneous change in the chemical indicator coloration from yellow to dark brown/black indicates correct air removal, adequate steam penetration and correct vacuum bomb operation. A failure can be evidenced as a lighter coloration, presence of spots or non-uniformity in the Test Sheet coloration.

## Indication for use

Chemdye® BDA4 Bowie-Dick Test Sheets are designed to be used in standard Bowie-Dick test packs assembled according to each country regulations. Test Sheets are used to test air removal efficiency at 132-134 °C in vacuum assisted steam sterilizers. The test pack should be placed horizontally on the lowest shelf in an already warmed, empty sterilizer. The sterilizer should run for 3.5 min. - 134 °C / 4 min. - 132 °C. Sterilizers showing non-uniform color change in the Chemdye® BDA4 Test Sheet should not be used and require technical service.

## Warning

Do not use the Chemdye® BDA4 Bowie-Dick Test Sheets for air removal monitoring of gravity air displacement steam sterilization processes, Dry Heat, Ethylene Oxide, or other low temperature sterilization processes.

## Storage

Store in a dark place and at temperatures between 10-30 °C, 30-80 % relative humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

It is recommended to protect the indicator Test Sheet after use inside a plastic foil for its storage.

## Directions for use

1. Perform the Bowie-Dick test daily and after sterilizer repair.
2. Prepare a test pack of absorbent surgical towels. The towels should be freshly laundered, but not ironed. The package should be conditioned in accordance with technical guidelines or in a way that its dimensions are at least 230 mm x 300 mm each side and 250-280 mm height.
3. Place a single Chemdye® BDA4 Bowie-Dick Test Sheet in the centre of the test pack.
4. By using self-adhesive tape, prepare a fabric standard package, without tightening excessively.
5. Place the test pack horizontally in the lower shelf of an empty vacuum assisted steam sterilizer.
6. Run the Bowie-Dick cycle for 3.5 min. - 134 °C / 4 min. - 132 °C.
7. After completion of the cycle, remove the Chemdye® BDA4 Bowie-Dick Test Sheet and examine:

## Satisfactory test

This is indicated by a Test Sheet that shows uniform dark brown/black color development. Such a test indicates rapid steam penetration and adequate air removal.

## Unsatisfactory test:

This is indicated by a heterogeneous color change into dark brown/black on the Test Sheet, such as lighter areas in the center. Unsatisfactory test should be reported immediately to the supervisor. The sterilizer should not be used until examined for possible malfunction.

Chemdye® BDA4 Bowie-Dick Test Sheet is manufactured according to ISO 11140 parts 1 and 3 standard.

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Terragene S.A