Multivariable Indicators

Result Reference Guide

Unprocessed


Failures

One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again.

The main failure causes are:

- Inadequate cycle duration
- Improper cycle selection
- Overloading / incorrect loading of the sterilizer
- Poor steam quality
- Sterilizer malfunction
- Blocking indicator with large metal mass

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NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

Explanation of Symbols

STEAM Product designed for use with steam sterilization cycles.

CE CE mark.
EC **REP** Authorized representative in the European Community.

MD Medical Device.

LOT Batch number.

EXP Expiration date.

M Manufacturer.

EN Multivariable Indicators

For steam sterilization processes

Type 4 indicators according to ISO 11140-1:2014

Product description

Chemdye® CD20 Multivariable Indicators have been designed to react to steam sterilization processes within loads, ensuring an adequate control of the effectiveness of sterilization processes (temperature, time and steam quality). Chemdye® CD20 Multivariable Indicators have been manufactured with an innovative ink technology, highly sensitive to steam quality and air presence within the load.

The purple ink was developed to turn to green when the process reaches the stated values of the critical process variables. Multivariable indicators allow to rapidly detect whether the critical parameters of the sterilization process for which they have been designed to respond have been reached.

Endpoint conditions: 121 °C, 15 minutes // 134 °C, 3.5 minutes.

Indicators 100 % free of toxic heavy metals.

Precautions

Do not expose Chemdye® CD20 Multivariable Indicators to ethylene oxide, dry heat or any sterilization processes other than steam.

Instructions for use

- Place a Chemdye® CD20 Multivariable Indicator within each package, bag or tray that will be sterilized by steam.
- In rigid containers place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.
- After completing the process, analyze results. The indicator must turn to the reference color for considering that the indicated conditions were met.

NOTE: if any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.

Shelf life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.

Endpoint stability reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard chemical indicators after use according to your country's healthcare and safety regulations.

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