## Multivariable Process Indicators

### Result Reference Guide

|  |  |
|--|--|
|  | <b>Unprocessed</b>   |
|  | <b>Failures</b>  |
|  | One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again. The main failure causes are:   |
|  | <ul style="list-style-type: none"> <li>• Inadequate cycle duration</li> <li>• Improper cycle selection</li> <li>• Overloading / incorrect loading of the sterilizer</li> <li>• Poor steam quality</li> <li>• Sterilizer malfunction</li> <li>• Blocking indicator with large metal mass</li> </ul> |
|  | <b>Processed and Satisfactory</b>  |
|  |  |

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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#### Explanation of Symbols

Product designed for use in Steam sterilization cycles

Batch number.

CE mark.

Expiration Date.

Authorized representative in the European Community.

Manufacturer.

### EN Multivariable Process Indicators

For Steam sterilization processes

Type 4 Indicators according to ISO 11140-1:2014

#### Product description

Chemdye® Multivariable Indicators have been designed to react within loads in Steam sterilization processes, ensuring an adequate control of the effectiveness of sterilization processes (temperature, time and steam quality). Blue reactive indicating ink has been developed to turn to black when the process reaches the stated values of the critical process variables. Endpoint conditions: 121 °C, 15 minutes // 134 °C, 3.5 minutes.

#### Characteristics

Punched strips of 200 x 15 mm, which allows to use the entire strip or half strip depending on the size of the package to be sterilized.

Indicators 100 % free of toxic heavy metals

#### Precautions

Do not expose Chemdye® Multivariable Indicators to any sterilization process other than Steam.

#### Instructions for use

• Place a Chemdye® Multivariable Indicator within each package, bag or tray that will be sterilized by Steam.

• In rigid containers place the indicator in every corner, or at least in two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.

• After completing the process analyze the results. The indicator must change to the reference color for considering that the indicated conditions were reached.

**NOTE:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

#### Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.

#### Shelf life

Indicators have an expiration date of 5 years from the date of manufacture when stored under the recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each product unit. Do not use indicators after their expiration date.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

#### Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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