## Multivariable Process Indicators

brown/black when the indicator is exposed to Steam considering the critical process variables which provides additional quality assurance. Endpoint condition: 121 °C, 15 minutes; 134 °C, 3.5 minutes.

Indicators 100 % free of toxic heavy metals.

### Precautions

Do not expose Chemdye® CD29-AD Multivariable Process Indicators to Ethylene Oxide, Dry Heat or any sterilization process other than Steam. Do not reuse.

### Instructions for Use

1. Place a Chemdye® CD29-AD Multivariable Process Indicator within each package, bag or tray that will be sterilized by Steam. In rigid containers, place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.

2. Perform the sterilization cycle.

3. After the sterilization process has finished, remove the indicators from the sterilizer and analyze the results. For all the chemical indicators exposed to Steam, the ink must turn to the dark brown/black reference color when the sterilization parameters are reached. For reference color, please refer to Result Reference Guide.

**NOTE:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

### Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store near sterilizing agents.

### Shelf Life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

### Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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

**Unprocessed**

### Failures


One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again.


The main failure causes are:

- Inadequate cycle duration
- Improper cycle selection
- Overloading / incorrect loading of the sterilizer

- Poor steam quality
- Sterilizer malfunction
- Blocking indicator with large metal mass

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**NOTE:** Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

### Explanation of Symbols

Product designed for use in Steam sterilization cycles.

Batch number.

CE mark.

Expiration Date.

Authorized representative in the European Community.

Manufacturer.

## EN Multivariable Process Indicators

For Steam sterilization processes

### Indications for Use

Chemdye® CD29-AD Chemical Indicator has been designed for monitoring Steam sterilization processes.

### Device Description

Chemdye® CD29-AD Multivariable Process Indicators (Type 4 according to ISO 11140-1:2014 standard) are single-use chemical indicators that consist of self-adhesive one-side laminated punched paper double strips printed with indicator ink. These indicators have been designed to monitor Steam sterilization processes within loads, ensuring an adequate exposure to the sterilizing agent during the sterilization process and allowing to distinguish between processed and unprocessed items.

Innovative highly sensitive yellow ink was developed to turn to dark

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