# CD40

Rev. 14 / 03.2024

## Multivariable Process Indicators

### Result Reference Guide

	<b>Unprocessed</b>
	<b>Failures</b> One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again. The main failure causes are: <ul style="list-style-type: none"> <li>Inadequate cycle duration</li> <li>Inadequate cycle selection</li> <li>Overloading / incorrect loading of the sterilizer</li> <li>Improper sterilization agent concentration</li> <li>Sterilizer malfunction</li> <li>Blocking indicator with large metal mass</li> </ul>
	<b>Processed and Satisfactory</b>
	<b>Processed and Satisfactory Long Exposure</b>

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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#### Explanation of Symbols

	Product designed for use with Plasma or Vaporized Hydrogen Peroxide sterilization cycles.		Medical Device.
	CE mark.		Batch number.
	Authorized representative in the European Community.		Expiration Date.
			Manufacturer.

#### Intended Use Table

Model	Vaporized Hydrogen Peroxide
	Cycles
CD40	STERRAD 100S 54 minutes
	STERRAD Standard and Express cycles
	V-Pro Max

Terragene® and Chemdye® are registered trademarks of Terragene S.A. STERRAD® is a registered trademark of Advanced Sterilization Products. V-PRO® is a registered trademark of Steris Corporation.

#### Trazanto

**ATTENTION:** If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto Automatic Quality Control and Traceability System** for cleaning, disinfection and sterilization processes. This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through **Bionova® Software**. Use Trazanto scanner to interpret the results (see instructions for use in the manual of the device). Do not write on the indicator before scanning.

## EN Multivariable Process Indicators

For Plasma or Vaporized Hydrogen Peroxide sterilization

#### Indications for Use United States

Terragene® Chemdye® CD40 is a Chemical Process Indicator intended for monitoring the efficacy of Vaporized Hydrogen Peroxide sterilization processes. The chemical indicator changes from purple to green or lighter to indicate that the conditions of the cycle have been met. **WARNING:** Contact with Hydrogen Peroxide during times and concentrations greater than those declared may cause the indicator to turn yellow. See Intended Use Table for reference.

#### Outside the United States

Chemdye® CD40 Chemical Indicator has been designed for monitoring Vaporized Hydrogen Peroxide sterilization processes.

#### Product Description

Chemdye® CD40 Multivariable Process Indicators (Type 4 according to ISO 11140-1:2014 standard) are single-use chemical indicators that consist of plastic strips printed with indicator ink. These indicators have been designed to monitor Plasma or Vaporized Hydrogen Peroxide sterilization processes within loads, ensuring an adequate exposure to the sterilizing agent during the sterilization process.

Innovative highly sensitive purple ink was developed to turn to green or lighter when the indicator is exposed to Plasma or Vaporized Hydrogen Peroxide considering the critical process variables which provides additional quality assurance.

Endpoint condition: 50 °C, 6 minutes, 2.3 mg/l H2O2.

**WARNING:** Contact with Hydrogen Peroxide during times and concentrations greater than those declared can cause the indicator ink to advance in the turn beyond its end point reaching yellow hues.

Indicators 100 % free of toxic heavy metals.

#### Precautions

Do not expose Chemdye® CD40 Multivariable Process Indicators to Ethylene Oxide, Dry Heat or any sterilization process other than Plasma or Vaporized Hydrogen Peroxide. Do not reuse.

#### Instructions for Use

1. Place a Chemdye® CD40 Multivariable Process Indicator within each package, bag or tray that will be sterilized by Plasma or Vaporized Hydrogen Peroxide. In rigid containers, place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.

2. Perform the sterilization cycle.

3. After the sterilization process has finished, remove the indicators from the sterilizer and analyze the results. For all the chemical indicators exposed to Plasma or Vaporized Hydrogen Peroxide, the ink must turn to the green or lighter reference color when the sterilization parameters are reached. For reference color, please refer to Result Reference Guide.

**NOTE:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

#### Monitoring frequency

Follow the policies and procedures with the monitoring frequency specified by the professional associations and/or standards corresponding to your country. As recommended practice, and to provide optimal patient safety, Terragene® recommends that each sterilization load be monitored with the appropriate biological indicator.

#### Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store near sterilizing agents.

#### Shelf Life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

#### Disposal

Discard chemical indicators after use along with plastic waste, according to your country's healthcare and safety regulations.

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