# CD42

Rev.1 / 05.2021

## Process Indicators

### Result Reference Guide

	Unprocessed
	Failures
	Processed and Satisfactory

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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#### Explanation of Symbols

	Product designed for use in Plasma or Vaporized Hydrogen Peroxide sterilization cycles.		Batch number.
	CE mark.		Manufacture Date.
	Authorized representative in the European Community.		Expiration Date.
			Manufacturer.

#### Intended Use Table

Model	Vaporized Hydrogen Peroxide
	Cycles
CD42	STERRAD 100S 54 minutes
	STERRAD Standard and Express cycles
	V-Pro Max and Sterizone VP4

Terragene® and Chemdye® are registered trademarks of Terragene S.A. STERRAD® is a registered trademark of Advanced Sterilization Products (a Johnson & Johnson Company). V-PRO® is a registered trademark of Steris Corporation. STERIZONE® is a registered trademark of TSO3 Inc

**ATTENTION:** If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto** Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes. This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through Bionova® Software. Use Trazanto scanner to interpret the results (see instructions for use in the manual of the device). Do not write on the indicator before scanning.

## Process Indicators

For Plasma or Vaporized Hydrogen Peroxide sterilization

### Indications for Use

**United States**  
Terragene® Chemdye® (CD42) is a Chemical Process Indicator intended for monitoring the efficacy of Vaporized Hydrogen Peroxide sterilization processes. The chemical indicator changes from red to yellow to indicate that the conditions of the cycle have been met. See Intended Use Table for reference.

### Outside the United States

Chemdye® CD42 Chemical Process Indicator has been designed for monitoring Vaporized Hydrogen Peroxide sterilization processes.

### Device Description

Chemdye® CD42 Process Indicators (Type 1 according to ISO 11140-1:2014 standard) are single-use chemical indicators that consist of plastic strips printed with indicator ink. These indicators have been designed to monitor Plasma or Vaporized Hydrogen Peroxide sterilization processes within loads, ensuring an adequate exposure to the sterilizing agent during the sterilization process and allowing to distinguish between processed and unprocessed items. Innovative highly sensitive red ink was developed to turn to yellow when the indicator is exposed to Plasma or Vaporized Hydrogen Peroxide. Endpoint condition: 50 °C, 6 minutes, 2.3 mg/l H2O2.

Indicators 100% free of toxic heavy metals.

### Precautions

Do not expose Chemdye® CD42 Process Indicators to Steam, Dry Heat or any sterilization process other than Plasma or Vaporized Hydrogen Peroxide. Do not reuse.

### Instructions for Use

- Place a Chemdye® CD42 Process Indicator within each package, bag or tray that will be sterilized by Plasma or Vaporized Hydrogen Peroxide. In rigid containers, place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.
- Perform the sterilization cycle.
- After the sterilization process has finished, remove the indicators from the sterilizer and analyze the results. For all the chemical indicators exposed to Plasma or Vaporized Hydrogen Peroxide, the ink must turn to the yellow reference color when the sterilization parameters are reached. For reference color, please refer to Result Reference Guide.

### Monitoring frequency

Follow facility Policies and Procedures which should specify a chemical indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a good practice and to provide optimal patient safety, Terragene® recommends that every sterilization load be monitored with an appropriate chemical indicator.

### Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % Relative Humidity. Do not wet. Do not store near sterilizing agents.

### Shelf Life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date. Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period 6 months when stored at previously indicated conditions.

### Disposal

Discard chemical indicators after use along with plastic waste, according to your country's healthcare and safety regulations.

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