
Chemical Indicators

Results Reference Guide


Unprocessed


Failures


Processed and Satisfactory

Processed and Satisfactory. Long exposure

Recommended Information

Sterilization date | Initials of the person in charge | Package content | Other information

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Explanation of Symbols

	Product designed for use in Plasma or Vaporized Hydrogen Peroxide sterilization cycles.		Medical Device.
	CE marked.		Batch number.
	Authorized representative in the European Community.		Expiration Date.
			Manufacturer.

EN Chemical Indicators

For Plasma or Vaporized Hydrogen Peroxide sterilization processes

Type 1 Indicators according to ISO 11140-1:2014

Product description

Chemdye® CD43 double self-adhesive labels have been designed to react to Plasma or Vaporized Hydrogen Peroxide sterilization processes within loads, ensuring an adequate exposure to the sterilizing agent during the sterilization process.

Chemdye® CD43 Indicators have been manufactured with an innovative highly sensitive ink, which was developed to turn from purple to green when the indicator comes into contact with the sterilizing agent. Its automatic record system allows control of the cycle by identifying the load, the sterilization date and, eventually, the sterilizer (in case of having more than one). Its self-adhesive properties offer the possibility of sticking it on the load or any other material to be sterilized. Chemdye® CD43 Indicators have been specially designed to adhere on non-woven wrappings, medical grade paper, polyester and polypropylene.

Characteristics

Double self-adhesive labels printed with indicator ink (29 x 22.5 mm), to be used with Chemdye® CG3 manual labeler.

Indicators 100 % free of toxic heavy metals.

Indications for use

Use Chemdye® CD43 indicators for monitoring Plasma or Vaporized Hydrogen Peroxide sterilization processes at 50 °C, 6 minutes, 2.3 mg/l H2O2. **WARNING:** Contact with Hydrogen Peroxide during times and concentrations greater than those declared can cause the indicator ink to advance in the turn beyond its end point reaching yellow hues.

Precautions

Do not use Chemdye® self-adhesive indicators to control Ethylene Oxide, Formaldehyde, Dry Heat or any other sterilization process different from Plasma or Vaporized Hydrogen Peroxide. Do not reuse chemical indicators.

Labeling using Chemdye® CG3 Labeler

The roll of Chemdye® CD43 indicators has a high-adhesion adhesive, specially designed for non-woven wrappings. If the application of the indicator label is performed using a CG3 labeler, pay attention to the following instructions:

1. Loading roll of labels in Chemdye® CG3 Labeler. Press the top black button (Fig. 1) and open the cover by lifting it completely until hearing a click. Do not force the cover once you have heard the click.
2. To load the labels roll, unroll the first 10 labels, detach and discard them, ensuring only the liner remains (Fig 2).
3. Place the labels roll between the fasteners as shown in (Fig. 3). The end of the roll should be at its top. Labels roll should be centered between the fasteners.
4. Pull the end of the roll and center it over the supply bar (Fig. 4). The material should emerge 15 to 16 cm.
5. Close the cover (Fig. 5). Be sure to hear the traction sound between the plastics when fully closed.
6. Press the trigger and insert the liner into the slot below the gray cylinder until it hits the top (Fig. 6).
7. The liner will loosen. To stretch it correctly, pull and release the trigger several times. The end of the liner should emerge from the back output of the labeler (Fig. 7).

B. Adjusting the print characters of Chemdye® CG3 Labeler.

The labeler offers the possibility of printing numbers in two lines, and letters and numbers in the remaining one. To make adjustments, each line has a knob.

1. Pull the knob out until the placeholder is on the desired position (Fig. 8).
2. Twist the knob until the desired character appears. It should be taken into consideration that for the two lines of numeric characters, the numbers to be printed are those seen in the setting window. For the alphanumeric line, attention should be drawn to the "Read or Print" conversion guide adhered to the labeler (Fig. 9), since the character to be printed is the opposite of the chosen one for this line. For example, in order to write "CG3" select, from the outside inwards, the character V (C will be printed), Z (G will be printed) and K (3 will be printed). In order to print a blank space, select the solid block for that position.
3. Repeat steps 1 and 2 for each change of position in the lines.
4. Push the knob in when you have finished setting the characters to print (Fig. 10).

C. Labeling of the load to be sterilized with Chemdye® CG3 Labeler.

Once loaded the roll and adjusted the printing characters, everything would be ready for labeling the load.

1. Press and release the trigger of the labeler as many times as needed until the label is exposed in the output position (Fig. 11). The trigger should be pressed hard and until it hits the top, thus preventing a defective feed of the liner and therefore an incorrect output of the label.
2. Place the labeler in vertical position to the labeling load (Fig. 12).
3. Put down the label on the labeling load. The printed side of the label must be visible and the adhesive side must contact the load (Fig. 12).
4. Slide the labeler in such a way that the label will be taped to the load until the entire label is exposed on the surface (Fig. 13).
5. The empty liner should emerge from the back output of the labeler (Fig. 14). This liner must be cut when its length hinders the labeling process.
6. Repeat steps 1 to 4 for each load to be labeled.

D. Changing the ink roller of Chemdye® CG3 Labeler.

Press the top black button (Fig. 15) and open the cover by lifting it completely until hearing a click. Do not force the cover once you have heard the click.

2. Unlock the fastener of the ink roller. Press down on the center portion of the fastener, where the arrow is (Fig. 16), until the ink roller is released from the fastener and can be removed.
3. In order to insert the new ink roller take it out of its wrap holding it at both ends in order to avoid staining. Place the ink roller in the labeler and press both ends firmly until it is fixed in the fastener (Fig. 17).
4. Close the cover (Fig. 18). Be sure to hear the traction sound between the plastics when fully closed.

E. Removal of the label roll from the Chemdye® CG3 Labeler.

1. Press the top black button (Fig. 19) and open the cover completely until you hear a click. Do not force the cover once you hear the click.
2. Remove the roll from the fasteners by moving the roll upward (Fig. 20).
3. Place the roll at the same level of the gray cylinder (Fig. 21).
4. Cut the tape in a location near to the gray cylinder (Fig. 21) and place the extracted roll in a safe place, away from the sterilizing agent.
5. Close the cover (Fig. 22). Be sure to hear the traction sound between the plastics when fully closed.
6. Remove the liner remnant from the slot on the back of the labeler by sliding it in the opposite direction of the labeler (Fig. 23).

F. Cleaning and clearing label jams in Chemdye® CG3 Labeler.

1. Press the top black button (Fig. 24) and open the cover completely until you hear a click. Do not force the cover once you hear the click.
2. Remove the roll from the fasteners by moving the roll upward (Fig. 25).
3. Place the roll at the same level of the gray cylinder (Fig. 26).
4. Cut the tape in a location near to the gray cylinder (Fig. 26) and place the extracted roll in a safe place, away from the sterilizing agent.
5. Remove the liner remnant from the slot on the back of the labeler by sliding it in the opposite direction of the labeler (Fig. 27).
6. Remove any loose labels on the inside of the labeler and under the supply bar (Fig. 28).
7. Use a commercial adhesive solvent (rather isopropyl alcohol) and a cotton swab to clean adhesive remnants.
8. Soak the swab in the solvent avoiding an excessive absorption.
9. Gently rub the swab onto surfaces and bearings where remaining adhesive is observed (Fig. 29).
10. Wait around 1.5 minutes the solvent evaporation.
11. Close the cover (Fig. 30). Be sure to hear the traction sound between the plastics when fully closed.

G. Clearing empty liner jam in the output slot of the Chemdye® CG3 Labeler.

If empty liner coming out is jammed during the labeling process:

1. Stop the labeling process.
2. Press the top black button (Fig. 31) and open the cover completely until you hear a click. Do not force the cover once you hear the click.
3. Remove the roll from the fasteners by moving the roll upward (Fig. 32).
4. Place the roll at the same level of the gray cylinder (Fig. 33).
5. Cut the tape in a location near to the gray cylinder (Fig. 33) and place the extracted roll in a safe place, away from the sterilizing agent.
6. Close the cover (Fig. 34). Be sure to hear the traction sound between the plastics when fully closed.
7. Press the labeler trigger while the liner remnant is removed in order to achieve a total clearance of the label (Fig. 35).
8. Visually check that no liner remains inside the labeler before it is reused. If remnants persist, shake the labeler (pressing and releasing the trigger) until the remnants falls (Fig. 36).

Sterilization

1. Sterilize as usual.
2. After the sterilization process has finished, open the sterilizer door, wait five minutes and remove the load. **CAUTION:** Wear safety glasses and gloves when removing the sterilized package.
3. Check the indicator ink printed on the label. A color change to green shows that the indicator was exposed to Plasma or Vaporized Hydrogen Peroxide. **IMPORTANT:** this color change does not fully evidence the process effectiveness to achieve sterility. If the chemical indicator color has not changed, check the sterilization process. **NOT:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Storage and shelf life

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents. The incorrect conservation of the chemical indicators may inactivate their functioning or alter the color after exposure to Steam.

Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each product roll. Do not use indicators after their expiration date. Chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard chemical indicators after use according to your country's healthcare and safety regulations.

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