CD47A

Rev. 0 / 10.2020

Chemical Indicator

EN Chemical Indicators

For monitoring environmental disinfection processes by Hydrogen Peroxide aerosol

Product description

Chemdye® CD47A chemical indicator has been specially developed to control environmental disinfection processes by Hydrogen Peroxide aerosol. It consists of a 15 mm diameter dot of self-adhesive synthetic substrate printed with reactive indicating ink. It undergoes a color change when the critical variables of the process reach specific values.

Instructions for use

1. Peel off Chemdye® CD47A chemical indicator from the roll and stick it on the surface which level of disinfection is pretended to be controlled. Use as many indicators as surfaces intended to be verified.
NOTE: it is recommended to stick Chemdye® CD47A chemical indicators on the surfaces that are more inaccessible for Hydrogen Peroxide aerosol, that are those ones that receive the lower flows, such as the corners of the room or the back of a bulky item or furniture.
2. Leave the room and run a normal disinfection cycle.
3. Once the process has finished, control the final color each Chemdye® CD47A chemical indicator used has reached. Each one must change from the initial purple color to green for the process to be considered successful.

Storage and shelf life

Store in a dry place, protected from light, at a temperature between 10-30 °C, and at a relative humidity between 30-80 %. Do not wet. Do not store close to disinfectant or sterilizing agents. Chemdye® CD47A chemical indicators have an expiration date of 5 years from the date of manufacture if stored under recommended conditions. Do not use Chemdye® CD47A chemical indicators after their expiration date.

Disposal

Discard used chemical indicators according to your country's healthcare and safety regulations.

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Terragene S.A.