CD83**Chemdye® ChemDose**

Rev. 3 / 11.2020

Chemical Indicator

EN**Chemical Indicator***For monitoring disinfection processes by continuous UV-C (254 nm) radiation*

It has been demonstrated that many microorganisms can be inactivated by continuous UV-C radiation. UV-C radiation is absorbed by DNA and RNA, causing changes in their structures and preventing their replication. The radiation dose required to achieve a microbial load reduction depends not only on the particular microorganism, but also on additional factors like radiation intensity (power) of the emitting source, distance between emitting source and surface, and time of exposure. Furthermore, surfaces that only receive radiation by reflection, get lower doses that may not be sufficient to inactivate the pathogenic microorganisms present. Considering all this, it is of utmost importance to control that all surfaces to be disinfected, receive the minimum radiation dose required.

Product Description

Chemdye® ChemDose CD83 chemical indicator has been specially developed to control disinfection processes by continuous UV-C radiation. It consists of a synthetic double self-adhesive substrate of 29 x 22.5 mm printed with reactive indicating ink. The reactive indicating ink has been calibrated to experience different color changes depending on the accumulated radiation dose received, allowing to quickly and accurately control the efficacy of the disinfection process.

Chemdye® ChemDose CD83 chemical indicators can be used either with or without Chemdye® CG3 Labeler. Using them with Chemdye® CG3 Labeler enables not only to speed-up the process of placing the indicators onto the surfaces which level of disinfection is pretended to be controlled but to record important information related to the disinfection cycle.

Instructions for use

When using Chemdye® ChemDose CD83 chemical indicator without Chemdye® CG3 Labeler

1. Peel off Chemdye® ChemDose CD83 chemical indicator from the roll, and stick it on the surface which level of disinfection is pretended to be controlled. Use as many indicators as surfaces intended to be verified.

NOTE 1: it is recommended to stick Chemdye® ChemDose chemical indicators onto plastic and fabric surfaces, in order to remove them easily after the disinfection process.

2. Run a normal disinfection cycle.

3. Once the process has finished, compare the final color reached by each Chemdye® ChemDose CD83 chemical indicator used against the

reference color printed on the label. If the final color matches the reference color or is redder than it, it can be considered that the surface has received at least a radiation dose of 50 mJ/cm².

NOTE 2: exposed and non-exposed Chemdye® ChemDose CD83 chemical indicators must be kept protected from light. Non-exposed Chemdye® ChemDose CD83 chemical indicators must be removed from their original packaging ONLY right before running the disinfection process.

When using Chemdye® ChemDose CD83 chemical indicator with Chemdye® CG3 Labeler

See Chemdye® CG3 Labeler Instructions for Use.

Storage and shelf life

Store in a dry place, protected from light, at a temperature between 10-30 °C, and at a relative humidity between 30-80 %. Do not wet. Do not store close to disinfectant or sterilizing agents.

Chemdye® ChemDose CD83 chemical indicators have an expiration date of 5 years from the date of manufacture if stored under recommended conditions.

Do not use Chemdye® ChemDose CD83 chemical indicators after their expiration date.

Disposal

Discard used chemical indicators according to your country's healthcare and safety regulations.

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└ Terragene S.A