# CD87-100

Rev. 7 / 02.2023

## Chemdye® CD87-100

### Chemical Indicator

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#### Explanation of Symbols

	CE mark.		Manufacture Date.
	Authorized representative in the European Community.		Expiration Date.
	Batch number.		Manufacturer.

### EN Chemical Indicator

For monitoring disinfection processes by continuous UV-C (254 nm) radiation

It has been demonstrated that many microorganisms can be inactivated by continuous UV-C radiation. UV-C radiation is absorbed by DNA and RNA, causing changes in their structures and preventing their replication. The radiation dose required to achieve a microbial load reduction depends not only on the particular microorganism, but also on additional factors like radiation intensity (power) of the emitting source, distance between emitting source and surface, and time of exposure. Furthermore, surfaces that only receive radiation by reflection, get lower doses that may not be sufficient to inactivate the pathogenic microorganisms present. Considering all this, it is of utmost importance to control that all surfaces to be disinfected, receive the minimum radiation dose required.

#### Product Description

Chemdye® ChemDose CD87-100 chemical indicator has been specially developed to control disinfection processes by continuous UV-C radiation. It consists of a 15 mm diameter dot of self-adhesive synthetic substrate printed with reactive indicating ink. The reactive indicating ink has been calibrated to experience different color changes depending on the accumulated radiation dose received, allowing to quickly and accurately control the efficacy of the disinfection process.

#### Instructions for use

1. Peel off Chemdye® ChemDose CD87-100 chemical indicator from the roll and stick it on the surface which level of disinfection is pretended to be controlled. Use as many indicators as surfaces intended to be verified.  
**NOTE 1:** It is recommended to distribute Chemdye® ChemDose CD87-100 chemical indicators uniformly along the room to be disinfected, placing at least one per square meter. It is recommended to place Chemdye® ChemDose CD87-100 chemical indicators not higher than 2 m over the floor level.
2. Leave the room and run a normal disinfection cycle.
3. Once the process has finished, contrast the final color of each Chemdye® ChemDose CD87-100 chemical indicator used against the Results Reference Card provided and determine the accumulated

radiation dose received. The attached Documentation Chart allows to quickly record and analyze the obtained results.

**NOTE 2:** Exposed and non-exposed Chemdye® ChemDose CD87-100 chemical indicators must be kept protected from light. Non-exposed Chemdye® ChemDose CD87-100 chemical indicators must be removed from its original packaging ONLY right before running the disinfection process.

**NOTE 3:** record the result immediately after the disinfection process.

**NOTE 4:** if any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established

#### Storage and shelf life

Store in a dry place, protected from light, at a temperature between 10-30 °C, and at a relative humidity between 30-80 %. Do not wet. Do not store close to disinfectant or sterilizing agents. Chemdye® ChemDose CD87-100 chemical indicators have an expiration date of 3 years from the date of manufacture if stored under recommended conditions.

Do not use Chemdye® ChemDose CD87-100 chemical indicators after their expiration date.

Once exposed, the indicator should be stored protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. The incorrect conservation of the chemical indicators may alter the color after exposure.

#### Disposal

Discard used chemical indicators according to your country's healthcare and safety regulations.

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