# CD87-1000M

Rev. 3 / 05.2021

 Chemdye® ChemDose

## Chemical Indicator

### Explanation of Symbols

	CE mark.		Manufacture Date.
	Authorized representative in the European Community.		Expiration Date.
	Batch number.		Manufacturer.

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### EN Chemical Indicator

For monitoring disinfection processes of N95 and FFP2 respirators by UVGI.

N95 and FFP2 filtering facepiece respirators (FFRs) are elements for respiratory protection which reuse has been studied in the context of the COVID-19 pandemic, by means of ultraviolet germicidal irradiation (UVGI) exposure. A report by the Center for Disease Control and Prevention (CDC) of the United States of America indicates that an UVGI dose of 1000 mJ/cm<sup>2</sup> can achieve a reduction of 99,9% of viruses similar to SARS-CoV-2. For this reason, Chemdye® ChemDose CD87-1000M chemical indicator has been calibrated to evaluate that FFRs have received the appropriate UVGI dose of 1000 mJ/cm<sup>2</sup>. According to the scientific consortium N95DECON, FFRs should not be reprocessed more than 5 times due to wear and tear. Chemdye® ChemDose CD87-1000M chemical indicator presentation allows to easily identify the amount of times a respirator has been reprocessed.

### Product Description

Chemdye® ChemDose CD87-1000M chemical indicator has been specially developed to control N95 y FFP2 respirators disinfection processes by UVGI. It consists of a 15 mm diameter dot of a self-adhesive synthetic substrate printed with reactive indicating ink. This ink has been calibrated to experience different color changes depending on the accumulated UVGI dose received, which allows to quickly and accurately control the efficacy of the disinfection process.

Chemdye® ChemDose CD87-1000M chemical indicator consists of five rolls numbered from 1 to 5, which allows the end-user to easily identify the amount of times their respirator has been reprocessed.

### Instructions for use

1. Before running the first disinfection process of the respirator, take off two Chemdye® ChemDose CD87-1000M chemical indicators from roll 1 and stick one of them onto the external surface and the other one onto the internal surface.
2. Run a normal disinfection process.
3. Once the process has finished, contrast the final color reached by both Chemdye® ChemDose CD87-1000M chemical indicators against the Reference Color Card provided in order to determine the UVGI dose received.

**NOTE 1:** in case one or none of the chemical indicators Chemdye® ChemDose CD87-1000M reaches the final color corresponding to the minimum UVGI dose required (1000 mJ/cm<sup>2</sup>), the disinfection process must be repeated following the indicated procedure in points 2.- to 3.- Otherwise, if both Chemdye® ChemDose CD87-1000M chemical indicators reach the final color corresponding to the minimum UVGI dose required, it can be considered the disinfection process has been effective and the respirator can continue to be used.

Before continue using the respirator, remove the Chemdye® ChemDose CD87-1000M chemical indicator that had been stucked onto the internal surface, but keep the one that had been stucked onto the external surface.

4. Before running the second disinfection process of the respirator, take off two Chemdye® ChemDose CD87-1000M chemical indicators from roll 2 and stick one of them onto the Chemdye® ChemDose CD87-1000M chemical indicator that had been used in the previous disinfection process and the other one onto the internal surface.

5. Run a normal disinfection process.

6. Once the process has finished, check the final color reached by both Chemdye® ChemDose CD87-1000M chemical indicators. A red final color corresponds to the minimum UVGI dose required of 1000 mJ/cm<sup>2</sup>.

**NOTE 2:** in case one or none of the chemical indicators Chemdye® ChemDose CD87-1000M reaches the final color corresponding to the minimum UVGI dose required (1000 mJ/cm<sup>2</sup>), the disinfection process must be repeated following the indicated procedure in points 5.- to 6.- Otherwise, if both Chemdye® ChemDose CD87-1000M chemical indicators reach the final color corresponding to the minimum UVGI dose required, it can be considered the disinfection process has been effective and the respirator can continue to be used.

Before continue using the respirator, remove the Chemdye® ChemDose CD87-1000M chemical indicator that had been stucked onto the internal surface, but keep the one that had been stucked onto the external surface.

7. Repeat the procedure indicated in points 4.- to 6.- for the third, fourth and fifth disinfection processes.

**NOTE 3:** exposed and non-exposed Chemdye® ChemDose CD87-1000M chemical indicators must be kept protected from light. Non-exposed Chemdye® ChemDose CD87-1000M chemical indicators must be removed from its original packaging ONLY right before running the disinfection process. Between two disinfection processes, the exposed Chemdye® ChemDose CD87-1000M chemical indicator that remains stucked onto the external surface of the respirator while it is being used, may change its color with respect to the one reached at the end of the disinfection process.

### Storage and shelf life

Store in a dry place, protected from light, at a temperature between 10-30 °C, and at a relative humidity between 30-80 %. Do not wet. Do not store close to disinfectant or sterilizing agents. Chemdye® ChemDose CD87-1000M chemical indicators have an expiration date of 5 years from the date of manufacture if stored under recommended conditions.

Do not use Chemdye® ChemDose CD87-1000M chemical indicators after their expiration date.

### Disposal

Discard used chemical indicators according to your country's healthcare and safety regulations.

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