CD87-200PX

Rev. 4 / 5.2021

Chemdy® ChemDose

Chemical Indicator

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Explanation of Symbols

CE mark.

EC REP Authorized representative in the European Community.

LOT Batch number.

Manufacture Date.

Expiration Date.

Manufacturer.

EN Chemical Indicator

For monitoring disinfection processes by pulsed UV (PX-UV) radiation

It has been demonstrated that many microorganisms can be inactivated by pulsed UV (PX-UV) radiation. Pulsed UV (PX-UV) radiation is absorbed by DNA and RNA, causing changes in their structures and preventing their replication. The radiation dose required to achieve a microbial load reduction depends not only on the particular microorganism, but also on additional factors like radiation intensity (power) of the emitting source, distance between emitting source and surface, and time of exposure. Furthermore, surfaces that only receive radiation by reflection, get lower doses that may not be sufficient to inactivate the pathogenic microorganisms present. Considering all this, it is of utmost importance to control that all surfaces to be disinfected, receive the minimum radiation dose required.

Product Description

Chemdy® ChemDose CD87-200PX chemical indicator has been specially developed to control disinfection processes by pulsed UV (PX-UV) radiation. It consists of a 15 mm diameter dot of self-adhesive synthetic substrate printed with reactive indicating ink. The reactive indicating ink has been calibrated to experience different color changes depending on the accumulated radiation dose received, allowing to quickly and accurately control the efficacy of the disinfection process.

Instructions for use

1. Peel off Chemdy® ChemDose CD87-200PX chemical indicator from the roll and stick it on the surface which level of disinfection is pretended to be controlled. Use as many indicators as surfaces intended to be verified.

NOTE 1: it is recommended to distribute Chemdy® ChemDose CD87-200PX chemical indicators uniformly along the room to be disinfected, placing at least one per square meter. It is recommended to place Chemdy® ChemDose CD87-200PX chemical indicators not higher than 2 m over the floor level.

2. Leave the room and run a normal disinfection cycle.

3. Once the process has finished, contrast the final color of each Chemdy® ChemDose CD87-200PX chemical indicator used against the

Results Reference Card provided and determine the accumulated radiation dose received.

NOTE 2: the Results Reference Guide has been calibrated against a particular pulsed UV radiation source, with a particular emission spectrum. There might be differences in the final colors reached at the different accumulated radiation doses received if a different pulsed UV radiation source is used. For more detailed information please contact Terragene®.

NOTE 3: exposed and non-exposed Chemdy® ChemDose CD87-200PX chemical indicators must be kept protected from light. Non-exposed Chemdy® ChemDose CD87-200PX chemical indicators must be removed from its original packaging ONLY right before running the disinfection process.

Storage and shelf life

Store in a dry place, protected from light, at a temperature between 10-30 °C, and at a relative humidity between 30-80 %. Do not wet.

Do not store close to disinfectant or sterilizing agents.

Chemdy® ChemDose CD87-200PX chemical indicators have an expiration date of 5 years from the date of manufacture if stored under recommended conditions.

Do not use Chemdy® ChemDose CD87-200PX chemical indicators after their expiration date.

Disposal

Discard used chemical indicators according to your country's healthcare and safety regulations.

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Terragene S.A