3D Chemical Indicators

Results Reference Guide

CDS47V	CDS47A
Unprocessed	
Failures	
Processed and Satisfactory	

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators. Images are reproduced for illustration purposes only, and do not necessarily match the article described in every detail.

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Explanation of Symbols

- MD** Medical Device.
- LOT** Batch number.
- ☞** Manufacture Date.
- 🕒** Expiration Date.
- 🏭** Manufacturer.
- CE** CE mark.
- EC REP** Authorized representative in the European Community.

EN Chemical Indicators
For monitoring of Airborne and Surface Disinfection processes with Hydrogen Peroxide Vapor or Aerosol

Product description
Chemical Indicators Chemdye® ChemSurf CDS47A and CDS47V have been designed to react in disinfection processes of airborne and surfaces by Hydrogen Peroxide. CDS47A and CDS47V inks were developed to change when the disinfection process reaches the values indicated for the critical process variables.

Characteristics

Chemdye® ChemSurf CDS47A and CDS47V Chemical Indicators consist of a pyramidal synthetic support (10.5 cm on the side and 9.0 cm in height) printed with a specific indicator ink designed to react to aerosol or vapor Hydrogen Peroxide, respectively. Chemdye® ChemSurf CDS47A and CDS47V Chemical Indicators present an exclusive design that allows it to dangle in height, according to the need. See forms of use in height and corresponding assemblies.

CDS47V (Vapor)

- Initial color: red
- Final color: yellow

CDS47A (Aerosol)

- Initial color: bordeaux
- Final color: green

Indicadores 100 % free of toxic heavy metals

Indications of use

Use Chemdye® ChemSurf CDS47A and CDS47V Chemical Indicators to monitor airborne and surface disinfection processes by vapor or aerosol Hydrogen Peroxide.

Conditions for CDS47V (Vapor)*

H₂O₂ solution concentration ≥ 35 %, temperature: 10-30 °C, relative humidity: 30-80 %.

Conditions for CDS47A (Aerosol)*

H₂O₂ solution concentration ≥ 6 %, temperature: 10-30 °C, relative humidity: 30-80 %.

(*) By individual enclosure.

Precautions

Do not use Chemdye® ChemSurf CDS47A and CDS47V Chemical Indicators to control disinfection processes with other agents than Hydrogen Peroxide or for sterilization processes by Hydrogen Peroxide in chamber.

Instructions for use

1. Perform the three-dimensional assembly of the indicators, following the assembly instructions specified below.
2. Place the indicators in those areas that you consider a priori more inaccessible to vapor or aerosol hydrogen peroxide. Locations that experience minimal H₂O₂ flows are suggested, such as the corners of the enclosure or behind a bulky item or furniture. To define which the most difficult areas to decontaminate are, it is recommended to use several indicators to validate the same disinfection cycle.
3. Perform the decontamination and airborne process as usual.
4. After the process is finished, analyze the results. The indicators must change to the reference color in all its faces to consider that the adequate disinfection conditions were reached.

NOTE: If the indicators are placed on a surface, the face in contact with it will not suffer a color change because it won't be exposed to the disinfecting agent.

IMPORTANT: The direct application of the sterilizing agent at a short distance, in a localized way and for a long time on the faces of the indicator can cause its damage.

IMPORTANT: Finger contact with reactive ink can alter the final color shift. Avoid excessive handling.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Storage and shelf life

Indicadores should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.

Indicadores have an expiration date of 5 years from the date of manufacture when stored under recommended conditions. Both the expiration date and lot number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date. Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard chemical indicators after use according to your country's healthcare and safety regulations.

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1. Fold the indicator by the dotted lines, leaving the chemical ink out.
2. Insert the flap "A" in the slot corresponding to the same letter, as indicated in Fig.1. Repeat the same for all the flaps.
3. To hang the indicator, pass the silicone seal through the perforations of faces 2 and 3. Fix the self-adhesive hook to a suitable surface and hang the pyramid from the seal.

Learn more about this product and check the digital color reference guide

Scan the QR code

CDS47A

CDS47V