CT22

Rev. 5/ 08/2022

Cintape®

## Process Indicators

### Result Reference Guide

Unprocessed

Failures

Processed and Satisfactory

**NOTE:** Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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#### Explanation of Symbols

 Product designed for use with Steam sterilization cycles.

 Batch number.

 CE mark.

 Manufacture Date.

 Authorized representative in the European Community.

 Expiration Date.

 Manufacturer.

#### Intended Use Table

Model	Autoclave/Steam Cycles		
	Gravity Displacement		
CT22	121 °C	132 °C	135 °C
	30 minutes	15 minutes 25 minutes	10 minutes
	Dynamic Air Removal (Vacuum Assist)		
	121 °C	132 °C	135 °C
	NA	4 minutes	3 minutes

Terragene® and Cintape® are registered trademarks of Terragene S.A.

## EN Process Indicators

For Steam sterilization

#### Indications for Use

##### United States

Terragene® Cintape® CT22 is a chemical process indicator tape intended for monitoring the efficacy of steam sterilization processes. The indicating tape changes from yellow to dark brown/black when exposure to steam. See Intended Use Table for reference.

##### Outside the United States

Cintape® CT22 is a chemical process indicator tape for external use for monitoring steam sterilization processes.

#### Device Description

Cintape® CT22 Process Indicators (Type 1 according ISO 11140-1:2014 standard) are single-use chemical indicators that consist of a roll of self-adhesive paper tape printed with indicator ink. These indicators have been designed to monitor Steam sterilization processes, ensuring an adequate exposure to the sterilizing agent during the sterilization process and allowing to distinguish between processed and unprocessed items.

The adhesive component of the tape allows the adhesion to different types of packaging and wraps, such as cloth, paper and plastic being easy to remove and leaving no residues. Innovative highly sensitive yellow ink was developed to turn to dark brown / black when the indicator is exposed to Steam. Endpoint condition: 121 °C, 10 minutes; 134 °C, 2 minutes.

#### Precautions

Do not expose Cintape® CT22 Process Indicators to Ethylene Oxide, Dry Heat or any sterilization process other than Steam. Do not reuse.

#### Instructions for Use

- Place within a sterilizer a wrapped package closed with Cintape® CT22 Process Indicator tape that will be sterilized by Steam.
- Perform the sterilization cycle.
- After the sterilization process has finished, remove the package from the sterilizer and analyze the indicator results. The indicator must turn to the reference color for considering that it was exposed to Steam. For reference color, please refer to Result Reference Guide.

**NOTE:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

#### Monitoring frequency

Follow facility Policies and Procedures which should specify a chemical indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a good practice and to provide optimal patient safety, Terragene® recommends that every sterilization load be monitored with an appropriate chemical indicator.

#### Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store near sterilizing agents.

#### Shelf Life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date. Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

#### Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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