Process Indicators

Result Reference Guide

	Unprocessed
	Failures
	Processed and Satisfactory

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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Explanation of Symbols

	Product designed for use with Plasma or Vaporized Hydrogen Peroxide sterilization cycles.		Medical Device.
	CE mark.		Batch number.
	Authorized representative in the European Community.		Manufacture Date.
			Expiration Date.
			Manufacturer.

Intended Use Systems

Model	Vaporized Hydrogen Peroxide
	Cycles
CT40	Sterrad® 100S Sterilization System
	Sterrad® NX Sterilization System (Standard and Advanced Cycles)
	Sterrad® NX with ALLClear® Technology Sterilization System (Standard and Advanced Cycles)
	Sterrad® 100NX Sterilization System (Standard, Flex, Express and Duo Cycles)
	Sterrad® 100NX with ALLClear® Technology Sterilization System (Standard, Flex, Express and Duo Cycles)
	V-Pro® maX Low Temperature Sterilization System (Non Lumen, Lumen and Flexible cycles)
	V-Pro® maX 2 Low Temperature Sterilization System (Fast Non Lumen, Non Lumen, Lumen and Flexible cycles)
	Sterizone® VP4 Sterilizer (Cycle 1)

Terragene® and Cintape® are registered trademarks of Terragene S.A.
STERRAD® is a registered trademark of Advanced Sterilization Products.
V-PRO® is a registered trademark of Steris Corporation.
STERIZONE® is a registered trademark of TSO₃ Inc.

EN Process Indicators

For Plasma or Vaporized Hydrogen Peroxide sterilization

Indications for Use

United States

Terragene® Cintape® Chemical Process Indicator Tape (CT40) is a chemical process indicator tape intended for monitoring the efficacy of vaporized hydrogen peroxide sterilization processes. The indicating tape changes from purple to green when exposure to vaporized hydrogen peroxide. See Intended Use systems for reference.

Outside the United States

Terragene® Cintape® CT40 is a chemical process indicator tape for external use for monitoring vaporized hydrogen peroxide sterilization processes.

Device Description

Terragene® Cintape® CT40 (Type 1 according ISO 11140-1:2014 standard) are single-use chemical indicators that consist of a roll of self-adhesive plastic tape printed with indicator ink. These indicators have been designed to monitor plasma or vaporized hydrogen peroxide sterilization processes, ensuring an adequate exposure to the sterilizing agent during the sterilization process and allowing to distinguish between processed and unprocessed items.

The adhesive component of the tape allows the adhesion to different types of packaging and wraps, such as cloth, paper and plastic being easy to remove and leaving no residues. Innovative highly sensitive purple ink was developed to turn to green when the indicator is exposed to plasma or vaporized hydrogen peroxide.

Indicators 100 % free of toxic heavy metals.

Precautions

Do not expose Cintape® CT40 to ethylene oxide (EO), dry heat or any sterilization process other than plasma or vaporized hydrogen peroxide. Do not reuse.

Instructions for Use

- Place within a sterilizer a wrapped package closed with Cintape® CT40 that will be sterilized by plasma or vaporized hydrogen peroxide.
- Perform the sterilization cycle.
- After the sterilization process has finished, remove the package from the sterilizer and analyze the indicator results. The indicator must turn to the reference color for considering that it was exposed to plasma or vaporized hydrogen peroxide. For reference color, please refer to the Result Reference Guide.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Monitoring frequency

Follow the policies and procedures with the monitoring frequency specified by the professional associations and/or standards corresponding to your country. As recommended practice, and to provide optimal patient safety, Terragene S.A. recommends that each sterilization lot be monitored with the appropriate biological indicator.

Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C (50-86 °F), 30-80 % relative humidity. Do not wet. Do not store near sterilizing agents.

Shelf Life

Chemical indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard chemical indicators after use along with plastic waste, according to your country's healthcare and safety regulations.

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