Unique point Chemical Integrator**Quality Certification****STEAM****Stated Values (SV)**

Temperature Time

LOT _____


121 °C _____ min.

W _____

128 °C _____ min.

X _____



135 °C _____ min.

SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

According to ISO 11140-1:2014 standard

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

Result Reference Guide**Unprocessed****Failures**

One or more parameters of sterilization were not achieved. The materials cannot be considered sterile and must be processed again. These results are usually obtained when steam quality is poor or adequate vacuum has not been generated.

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Explanation of Symbols

STEAM Product designed for use with Steam sterilization cycles.

MD Medical Device.

CE mark.

LOT Batch number.

EC REP Authorized representative in the European Community.

Expiration Date.

Manufacturer.

Intended Use Table

| Model | Autoclave/Steam Cycles | | |
|--|------------------------|-------------|------------|
| | Gravity Displacement | | |
| 121 °C | 132 °C | 135 °C | |
| 30 minutes | 15 minutes | 10 minutes | 25 minutes |
| Dynamic Air Removal (Vacuum Assist) | | | |
| IT26-1YS | 121 °C | 132 °C | 135 °C |
| NA | 4 minutes | 3 minutes | |
| Minimum Stated Values (SV) for Chemical Integrators | | | |
| 121 °C | 132 °C | 135 °C | |
| 16.5 minutes | 2.0 minutes | 1.2 minutes | |

Terragene® and Integron® are registered trademarks of Terragene S.A.

ATTENTION: If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive Trazanto Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes.

This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through Bonivona® Software. Use Trazanto scanner to interpret the results (see instructions for use in the manual of the device).

Do not write on the indicator before scanning.

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