Chemical Integrator



Rev. 15 / 04.2024




IT26-C | IT26-C Extender

Result reference guide


Unprocessed

Accept

Reject

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Uso exclusivo para profesionales e Instituciones Sanitarias
Technical Director / Director Técnico. Lic. Adrián J. Rovetto

NOTE: The extender of IT26-C Extender version is intended to help retrieving the indicator from the pouch/packs.

Explanation of Symbols

STEAM Product designed for use with Steam sterilization cycles.

MD Medical Device.

LOT Batch number.

CE mark.

M Manufacture Date.

Authorized representative in the European Community.

E Expiration Date.

M Manufacturer.

Intended Use Table

Model	Autoclave/Steam Cycles		
	Gravity Displacement		
IT26-C	121 °C 30 minutes	132 °C 15 minutes 25 minutes	135 °C 10 minutes
	Dynamic Air Removal (Vacuum Assist)		
	121 °C NA	132 °C 4 minutes	135 °C 3 minutes
	Minimum Stated Values (SV) for Chemical Integrators		
	121 °C 16.5 minutes	132 °C 2.0 minutes	135 °C 1.2 minutes

Trazanto

ATTENTION: If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto** Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes.

This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through **Bionova® Software**. Use **Trazanto** scanner to interpret the results (see instructions for use in the manual of the device).

ATTENTION: do not write on the indicator before scanning.

Terragene® and Integron® are registered trademarks of Terragene S.A.
ATCC® is a registered trademark of American Type Culture Collection.

EN Chemical Integrator

For Steam sterilization

Indications for Use

United States

The integrator Terragene® Integron® IT26-C is designed to chemically react over time with the critical parameters of steam sterilization cycle within a specified tolerance. The integrating indicator strip is intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for the following sterilization cycles: Gravity Displacement: 121 °C for 30 minutes, 132 °C for 15 minutes, 132 °C for 25 minutes, 135 °C for 10 minutes; Dynamic Air Removal (Vacuum Assist): 132 °C for 4 minutes, 135 °C for 3 minutes. SV121 °C/16.5 min, SV132 °C/2.0 min, SV135 °C/1.2 min. See Intended Use Table for reference.

Outside the United States

Integron® IT26-C moving front chemical integrator is designed for monitoring Steam sterilization processes between 118 °C (244.4 °F) and 138 °C (280.4 °F).

Device Description

Integron® IT26-C Chemical Integrator (Type 5 according to ISO 11140-1:2014 standard) are single-use moving front chemical indicators for monitoring steam sterilization processes. These indicators ensure an adequate control of the effectiveness of steam sterilization processes by monitoring all critical parameters. Integron® IT26-C integrator indicator consists of a special paper wick and a steam and temperature sensitive chemical pellet contained in a paper/film/foil laminate. The chemical pellet melts and migrates as a dark bar along the paper wick. The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met. The extent of migration depends on steam quality, time and temperature. The ACCEPT result is reached when a theoretical spore population reaches its kill time, indicating integration condition

has been reached. This condition is calibrated with the kill time of a 10⁶ *Geobacillus stearothermophilus* ATCC® 7953 spores population, calculated in BIER (Biological Indicator Evaluator Resistometer).
NOTE: The extender of IT26-C Extender version is intended to help retrieving the indicator from the pouch/packs.

Precautions

Do not expose Integron® IT26-C to EO, Dry Heat, or any sterilization process other than Steam. Do not reuse.

Instructions for Use

1. Carefully open the package and remove only the number of indicators needed. Reseal the package tightly.
 2. Place an Integron® IT26-C indicator in each pack, peel pouch, container system or tray to be steam sterilized in the most inaccessible area to sterilizing agent.
 3. Perform of sterilization cycle.
 4. After the sterilization process has finished, the dark bar should have reached the ACCEPT zone. If the dark bar has not reached the ACCEPT zone, a REJECT result is indicated and the items in the pack, peel pouch, container system, or tray were not exposed to proper steam sterilization conditions. These items should be reprocessed according to established procedures and a new IT26-C has to be used.
- NOTE: If a damp appearance is observed on the IT26-C indicator after the sterilization cycle, Terragene® recommends repeating the sterilization cycle. If the problem persists, Terragene® recommends scheduling equipment maintenance.
NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Monitoring frequency

Follow facility Policies and Procedures which should specify a chemical indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a good practice and to provide optimal patient safety, Terragene® recommends that every sterilization load be monitored with an appropriate chemical indicator.

Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C (50-86 °F), 30-80 % relative humidity. Do not get wet. Do not store near sterilizing agents.

Shelf Life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Both the expiration date and batch number of the product are specified on the package label. Do not use indicators after the expiration date.
Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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