Sequential integrator indicator

Quality Certification

STEAM

Product code _____


LOT _____

Temperature _____ Time (min.) _____

121 °C _____

128 °C _____

135 °C _____

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Type 5 Indicator according to ISO 11140-1:2014

Product description

Integron® IT26 are three-level sequential integrating indicators designed to react inside loads exposed to Steam sterilization processes. The Integration level consists of a pink ink developed to turn to black when a theoretical population of spores reaches its kill time, indicating that the **SAFE** condition is achieved.

Level 1 - **Blue: Exposure level**, it turns to black to indicate the exposure to the process but it does not assure that the conditions indicated in the integrator have been met (**UNSAFE** condition).

Level 2 - **Pink: Integration level**, it turns to black thus indicating that conditions indicated in the integrator have been met. These conditions are calibrated with the kill time of a 10⁶ theoretical *Geobacillus stearothermophilus* spore population, calculated in a BIER (Biological Indicator Evaluator Resistometer). This is called **SAFE** condition.

Level 3 - **Green: Over-integration level**, it turns to black when integration conditions indicated in the product have been overreached. This is called **EXTENDED** condition.

Characteristics

The indicator consists of a one-side-laminated paper strip of 60 mm x 25 mm.

Indications for use

Use integrator indicators for monitoring Steam sterilization processes. **SAFE** condition is reached in Steam sterilization processes as described in the quality certificate.

Precautions

Do not expose this product to EO, Dry Heat, or any sterilization processes other than Steam. Do not reuse the sterilizer until the integrator's color reaches or passes the **SAFE** condition.

Instructions of use

- Place an Integron® IT26 Integrator Indicator within each package, bag or tray that will be sterilized by Steam.
- In rigid containers place an Integron® IT26 Integrator Indicator in every corner, or at least two diagonally or opposite corners, to ensure that steam has reached the most inaccessible areas.
- After completing the process, analyze the results. If the **SAFE** integration condition has been reached or overreached, the sterilization process was successful.

Storage and shelf life

- Indicators should be stored in dry places, away from light, at temperatures between 10-30 °C, 30-80 % Relative Humidity. Do not wet. Do not store close to sterilizing agents.

- Indicators have an expiration date of 5 years from the date of manufacture when stored in the recommended conditions. Both expiration date as well as the product batch number are specified on the package label and on each unit of product. Do not use indicators after their expiration date.

- After use, the indicator will not change in 6 months if stored in temperature and relative humidity previously indicated.

Disposal

Discard along with cardboard waste according to your country's healthcare regulations.

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