Unique point Integrators

Rev. 12 / 04.2024

Integron®

IT26S | IT26SBL | IT26SAD

Quality Certification


STEAM

Product code _____

LOT

Stated Values (SV)

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SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

According to ISO 11140-1:2014 standard

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

Result reference guide


Unprocessed

Failures

One or more parameters of sterilization were not achieved. The materials cannot be considered sterile and must be processed again. These results are usually obtained when steam quality is poor or adequate vacuum has not been generated.

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Explanation of Symbols

Product designed for use with Steam sterilization cycles.

Medical Device.

Batch number.

CE mark.

Expiration Date.

Authorized representative in the European Community.

Manufacturer.

EN Unique point Integrator Indicators

For Steam sterilization

Type 5 Indicators according to ISO 11140-1:2014

Product description

Integron® unique point Integrator Indicators were developed for monitoring Steam sterilization processes between 121 °C and 135 °C, ensuring an adequate control of the effectiveness of sterilization processes (temperature, time, quality of steam). Purple ink was developed to turn to green when a theoretical spores population reaches its kill time, indicating that integration condition has been reached.

This condition is calibrated with the kill time of a 10⁶ *Geobacillus stearothermophilus* spore population, calculated in a BIER (Biological Indicator Evaluator Resistometer).

Integrator Indicators are manufactured with inks that are highly sensitive to the Steam sterilization process variables.

Characteristics

IT26S: paper strip of 72 mm x 23 mm.

IT26SBL: two-sided laminated paper strip of 72 mm x 23 mm.

IT26SAD: paper strip of 72 mm x 23 mm in a self-adhesive support for record.

Indicators 100 % free of toxic heavy metals.

Indications for use

Use the integrator indicator to monitor Steam sterilization processes between 121 °C and 135 °C. Integration condition is reached in Steam sterilization processes as described in the quality certificate.

Precautions

Do not expose this product to Ethylene Oxide, Dry Heat, or any chemical or physical sterilization processes other than Steam.

Instructions of use

1. Place an Integron® unique point Integrator Indicator within each package, bag or tray that will be sterilized by Steam.

2. In rigid containers place an Integron® unique point integrator indicator in every corner, or at least in two diagonally opposite corners, to ensure that steam has reached the most inaccessible areas.

3. After completing the process, analyze the results. If the ink printed on the indicator changes to the reference color printed on the product, integration condition is met and the sterilization process was successful.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Storage

Indicators should be stored in dry places, away from light, at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store near sterilizing agents.

Shelf life

Indicators have an expiration date of 5 years from the date of manufacture when stored under recommended conditions. The expiration date as well as the batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date. The end point of the chemical indicator shall remain unchanged for a period of not less than 6 months when stored under the conditions specified above.

Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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