Unique point Emulator Indicator

IT27-3.5S | IT27-3.5SAD | IT27-5.3S | IT27-5.3SAD
IT27-7S | IT27-7SAD | IT27-9S

Result Reference Guide

Unprocessed

Failures


One or more parameters of sterilization were not achieved. The materials cannot be considered sterile and must be processed again. These results are usually obtained when steam quality is poor or adequate vacuum has not been generated

Processed and satisfactory

NOTE: Reference colors exhibited in printed prospectus as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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EN Unique point Emulator Indicators

For monitoring steam sterilization processes

Type 6 indicator according to ISO 11140-1:2014

Product description

Integron® Emulator Indicators have been designed to react in steam sterilization processes when placed inside of the loads, ensuring an adequate control of the effectiveness of sterilization processes (temperature, time, steam quality). The purple reactive indicating ink has been developed to change to green when the process reaches the established values of the critical parameters.

Characteristics

Integron® IT27-3.5S, IT27-3.5SAD, IT27-5.3S, IT27-5.3SAD, IT27-7S, IT27-7SAD and IT27-9S Emulator Indicators consist of a paper strip of 70 mm x 23 mm. Integron® IT27-3.5SAD, IT27-5.3SAD and IT27-7SAD Emulator Indicators are printed onto a self-adhesive support which allows their easy record.

Indicators are 100% free of toxic heavy metals.

Indications for use

Use Integron® Emulator Indicators for monitoring steam sterilization processes.

Emulsion condition:

IT27-3.5S	3.5 min. 134 °C	
IT27-3.5SAD	3.5 min. 134 °C	
IT27-5.3S	5.3 min. 134 °C	15 min. 121 °C
IT27-5.3SAD	5.3 min. 134 °C	15 min. 121 °C
IT27-7S	7 min. 134 °C	20 min. 121 °C
IT27-7SAD	7 min. 134 °C	20 min. 121 °C
IT27-9S	9 min. 134 °C	

Precautions

Do not expose these products to sterilization processes other than Steam. Do not reuse the sterilizer until Integron® Emulator Indicators reach the final color.

Instructions of use

- Place an Integron® Emulator Indicator within each package, bag or tray that will be sterilized by Steam.
- In rigid containers place an Integron® Emulator Indicator in every corner, or at least two diagonally or opposite corners, to ensure the steam will reach the most inaccessible areas.
- After completing the process, analyze the results. If Integron® Emulator Indicator changes to the reference color it can be considered that emulsion conditions have been reached.

Storage

Indicators should be stored in a dry place, protected from light, at a temperature between 10-30 °C, and at a relative humidity between 30-80 %. Do not wet. Do not store close to sterilizing agents.

Shelf life

Integron® Emulator Indicators have an expiration date of 5 years from the date of manufacture when stored under recommended conditions. Both expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use Integron® Emulator Indicators after their expiration date.

Endpoint Stability Reaction: Integron® Emulator Indicators endpoint shall remain unchanged for a period of 6 months when stored at previously indicated conditions.

Disposal

Discard along with cardboard waste according to your country's healthcare regulations.

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