Emulator Indicator

Results Reference Guide

Unprocessed

Failures

One or more parameters of the sterilization cycle were not achieved. The materials cannot be considered sterile and must be processed again.

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Explanation of symbols

STEAM

Product designed for use with steam sterilization cycles.

MD

Medical Device.

CE

CE mark.

L01

Batch number.

EC REP

Authorized representative in the European Community.

X

Expiration date.

M

Manufacturer.

EN Emulator Indicator

For steam sterilization

Type 6 Indicator according to ISO 11140-1:2014

Product Description

Integron® Emulating Indicators have been designed to react to steam sterilization processes inside loads. The blue ink was designed to turn to black when the process reaches the stated values (SV) of the critical

process variables indicated in the emulator.

Characteristics

The IT28-AD indicator consists of 60 mm x 25 mm single sided laminated self-adhesive paper printed.

Indicators 100 % free of toxic heavy metals

Indications for Use

Use the emulator indicator for monitoring steam sterilization processes. The endpoint is reached in steam sterilization processes of 3.5 minutes at 134 °C or 15 minutes at 121 °C.

Precaution

Do not use this emulator indicator for monitoring ethylene oxide, dry heat or other sterilization processes different to steam.

Instructions of Use

1. Place an Integron® Indicator Emulator within each package, pouch or tray to be sterilized by steam.
2. In rigid containers place an Integron® Emulator Indicator in every corner, or at least two diagonally or opposite corners, to ensure that the steam has reached the most inaccessible areas.
3. After completing the process, analyze the results. If the emulsion condition has been reached, the indicator ink must turn to black and it is concluded the sterilization process was successful.

NOTE: if any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Storage

Indicators should be stored in dry places, away from light, at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.

Shelf Life

Indicators have an expiration date of 5 years from the date of manufacture when stored in the recommended conditions. Both expiration date as well as the batch of the product are specified on the package label and on each unit of product. Do not use indicators after their expiration date.

Endpoint stability reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard along with cardboard waste according to your country's healthcare regulations.

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