# IT31

Rev. 4 / 07.2021

**Integron**

## Emulator Indicator

### Results Reference Guide

#### Unprocessed

#### Failures

One or more parameters of sterilization were not achieved. The materials cannot be considered sterile and must be processed again.

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#### Explanation of Symbols

**DRY** Product designed for use in Dry Heat sterilization processes.

**LOT** Batch number.

**Exp** Expiration Date.

**M** Manufacturer.

**CE** CE mark.

**EC REP** Authorized representative in the European Community.

**CE**

## EN Emulator Indicator

For Dry Heat sterilization processes

Type 6 Indicator according to ISO 11140-1:2014

#### Product description

Integron® IT31 Emulator Indicator has been designed to react to Dry Heat sterilization processes at 160 °C, 40 minutes. Its blue indicating ink has been developed to turn to brown when the process reaches the stated value (SV) of the critical process variables for which it has been designed.

#### Characteristics

IT31 Emulator Indicator consists of a paper strip of 60 mm x 25 mm printed with indicating ink.

Indicadores 100 % free of toxic heavy metals

#### Indications for use

Use Integron® IT31 Emulator Indicators to monitor Dry Heat sterilization processes at 160 °C, 40 minutes.

#### Precautions

Do not expose this product to EO, Steam, or any chemical or physical agent sterilization processes other than Dry Heat. Do not reuse the sterilizer until the emulator reaches its endpoint.

#### Instructions of use

1. Place an emulator indicator within each package or tray that will be sterilized by Dry Heat.

2. After completing the process, analyze results.

**NOTE:** if any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

#### Interpretation of the results

If the indicator has reached its endpoint, the sterilization process was successful.

Initial color: blue

Final color: brown

#### Storage

Indicators should be stored in dry places, away from light at temperatures between 10-30 °C, 30-80 % relative humidity.

Do not wet.

Do not store close to sterilizing agents.

Do not reuse.

#### Shelf life

Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. The expiration date as well as the batch number are specified on the package label and in each product unit.

Do not use indicators after their expiration date.

Chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

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