
# Helix-PCD

Rev. 8/ 07.2024

## Result reference guide

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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### Explanation of Symbols

	Product designed for use in Ethylene Oxide sterilization processes.		Medical Device.
	CE mark.		Batch number.
	Authorized representative in the European Community.		Expiration Date.
			Manufacturer.

## EN PCD-A-E1

For monitoring ethylene oxide sterilization processes

Chemical strips (Type 2 Chemical Indicator according to ISO 11140-1:2014) for use with Helix-PCD in the control of ethylene oxide sterilization processes.

Indicator 100 % free of toxic heavy metals

### Instructions for use

- Insert PCD-A-E1 chemical indicator, according to the selected process conditions, within the pin (metal capsule) or slot (plastic capsule) of the screw cap; settle the reactive indicating ink towards the capsule lumen, avoiding folding those areas printed with reactive indicating ink.
- IMPORTANT:** avoid contact between PCD-A-E1 chemical indicator and capsule interior surface, because it can alter the result.
- Put the capsule cap with moderate pressure.
- Settle the system into the fabric bag.
- The entire system should be sterilized according to the usual cycle used. For example: 240 minutes, 600 mg/l, 54 °C, 60% RH. Do not use cycles below 90 minutes and EO 300 mg/l.  
Initial color: sky blue  
Final color: green
- After the period of ventilation, remove the indicator and evaluate the result. If the 4 bars change to the final color, the result can be considered successful, which means the sterilizer is working correctly. If one bar doesn't reach the final color, it means that residual air is present in the sterilizer, therefore a malfunction of the equipment. We suggest to perform the test once again, before sending the equipment to technical service.

NOTE: if any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

### Storage

Store in a dry place, protected from light, at a temperature between 10-30 °C, and at a relative humidity between 30-80 %. Do not wet. Do not store close to sterilizing agents.

### Shelf life

PCD-A-E1 chemical indicators have an expiration date of 5 years from the date of manufacture when stored under recommended storage conditions. Both, expiration date and batch number of the product, are specified on the package label and in each product unit.  
Do not use PCD-A-E1 chemical indicators after their expiration date.  
Do not reuse.  
Endpoint reaction stability: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

### Disposal

Discard used chemical indicators along with paper waste, according to your country's healthcare and safety regulations.

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