## Result Reference Guide


**NOTE:** Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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#### Explanation of Symbols

 VH202	Product designed for Hydrogen Peroxide sterilization processes.	 MD	Medical Device.
 CE	CE mark.	 LOT	Batch number.
 EC REP	Authorized representative in the European Community.	 Exp. Date.	Expiration Date.
		 Manufacturer.	Manufacturer.

### EN Helix-PCD

For Hydrogen Peroxide sterilization

Chemical Strips (Type 2 Chemical Indicator according to ISO 11140-1) for use in Helix-PCD for Hydrogen Peroxide.

#### Instructions for use

1-A. Remove the threaded cap from the metal or plastic capsule. Only for metal capsule Helix devices, remove the bolt attached to the cap (the bolt should NOT be used in Hydrogen Peroxide processes).

1-B. Place of chemical indicator in the capsule. Settle the reactive ink toward the capsule lumen, avoiding folds into the chemical ink. **IMPORTANT:** Avoid the contact between the chemical indicator and the capsule internal surface, since the result can be altered.

2. Place the lid of the capsule with moderate pressure. **ATTENTION:** Do not place the system inside the cloth bag, since it can interfere with the process. The cloth bag must be used in steam processes.

**NOTE:** Terragene® suggests using the system Chemdye® Helix-PCD designed for small sterilizers (CDPCD2X025).

3. Start the selected program and test the process.

Initial color: purple

Final color: green

4. After the period of ventilation, remove the indicator and evaluate the result. If the 4 bars changed to the final color, the result is successful, the sterilizer is working correctly. If one bar is different to the final color, it shows the presence of residual air in the sterilizer, therefore an equipment malfunction.

We suggest retrying the test once again, before sending the equipment to technical service.

**NOTE:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

#### Storage and shelf life

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.

Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.

Endpoint stability reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

#### Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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