Helix-PCD


Rev. 5 / 06.2020

Result reference guide

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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EN PCD-A-F1

For Monitoring of Formaldehyde Sterilization Processes

Chemical Strips (Type 2 Chemical Indicator according to ISO 11140-1:2014) for use with Helix-PCD in the monitoring of Formaldehyde Sterilization Processes.

Instructions for use

1. Insert the chemical indicator according to the selected process conditions within the pin (metal capsule) or slot (plastic capsule) of the screw cap; settle the reactive indicating ink towards the capsule lumen, and avoid folding zones printed with reactive indicating ink. **IMPORTANT:** Avoid contact between chemical indicator and capsule interior surface, because result can be altered.
2. Put the capsule cap with a moderate pressure. **CAUTION!** Do not put the system into the fabric bag, it could interfere in the process.
3. For LTSF processes, start the selected program and test the process.
Initial color: bordeaux
Final color: green
4. After the period of ventilation, remove the indicator and evaluate the result. If the 4 bars reach the final color, the result can be considered successful, the sterilizer is working correctly. If one bar doesn't reach the final color, it means that residual air is present in the sterilizer, therefore a malfunction of the equipment. We suggest running the test once again, before sending the equipment to technical service.

Storage

Indicators should be stored in dry places, protected from light, at temperatures between 10-30 °C, and at a relative humidity between 30-80 %.
Do not wet.
Do not store close to sterilizing agents.
Do not use indicators after their expiration date.
Do not reuse.

Shelf life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both the expiration

date and batch number of the product are specified on the package label and in each unit of the product.
Endpoint reaction stability: chemical indicator endpoint shall remain unchanged for a period of 6 months when stored at previously indicated conditions.

Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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Terragene S.A