
PCD-A-3.5Y | PCD-A-5.3Y | PCD-A-7.0Y
PCD-A-9.0Y | PCD-A-18.0Y

Result Reference Guide


Temperatura, tiempo y penetración de vapor suficientes / Sufficient temperature, time and steam penetration / Temperatura, tempo e penetração do vapor suficientes / Temperature, temps et pénétration de vapeur suffisants / Sufficiente temperatura, il tempo e la penetrazione del vapore / Yeterli sıcaklık, zaman ve buhar nüfuzu.

Remoción de aire y penetración de vapor insuficientes / Insufficient air removal and steam penetration / Remoção do ar e penetração do vapor insuficientes / Enlèvement d'air et pénétration de vapeur insuffisants / Insufficiente rimozione dell'aria e penetrazione del vapore / Yetersiz hava giderilmesi ve buhar nüfuzu.

Temperatura alcanzada. Sin remoción de aire ni penetración de vapor / Temperature achieved, but no air removal and no steam penetration / Temperatura se alcanza. Sem remoção do ar ou penetração de vapor / Temperature obtenue. Sans enlèvement d'air ni pénétration de la vapeur / Temperatura raggiunta, ma non la rimozione dell'aria e la penetrazione del vapore / Sıcaklığa ulaşılmış, ancak hava giderilmesi ve buhar nüfuzu yok.

Temperatura insuficiente. Sin remoción de aire ni penetración de vapor / Insufficient temperature and no air removal and no steam penetration / Temperatura insuficiente. Sem remoção do ar ou penetração de vapor / Temperature insuffisante. Sans enlèvement d'air ni pénétration de la vapeur / Temperatura insufficiente e non rimozione dell'aria e penetrazione del vapore / Yetersiz sıcaklık ve hava giderimi yok ve buhar nüfuzu yok.

NOTE: Reference colors exhibited in printed prospects and boxes as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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Explanation of Symbols

Product designed for use with Steam sterilization cycles.

CE mark.

Authorized representative in the European Community.

Batch number.

Expiration Date.

Manufacturer.

EN Helix-PCD

For monitoring of steam sterilization processes

Chemical strips (Type 2 Chemical Indicator according to ISO 11140-1:2014) for use with Helix-PCD in the control of steam sterilization processes.

Indicators are 100% free of toxic heavy metals.

Instructions for use

1. Insert PCD-A chemical indicator, according to the selected process conditions, within the holder (metal capsule) or slot (plastic capsule) of the screw cap; settle the reactive indicating ink towards the capsule lumen, avoiding folding those areas printed with reactive indicating ink.

IMPORTANT: avoid contact between PCD-A chemical indicator and capsule interior surface, because it can alter the result.

2. Put the capsule cap with moderate pressure.

3. Settle the system into the fabric bag.

4. The entire system should be sterilized in a vacuum chamber under conditions that depend on the chemical indicator to be used:

PCD-A-3.5Y 3.5 min. 134 °C

PCD-A-5.3Y 5.3 min. 134 °C, 15 min. 121 °C

PCD-A-7.0Y 7.0 min. 134 °C, 20 min. 121 °C

PCD-A-9.0Y 9.0 min. 134 °C

PCD-A-18.0Y 18.0 min. 134 °C

5. Remove the indicator and evaluate the result. If the 4 bars reach the final color, the result can be considered successful, which means the sterilizer is working correctly. If one bar doesn't reach the final color, it means that residual air is present in the sterilizer, therefore a malfunction of the equipment. The reasons for the failure can be insufficient air removal, leakages and noncondensable gases in steam. We suggest to perform the test once again, before sending the equipment to technical service.

Storage

Store in a dry place, protected from light, at a temperature between 10-30 °C, and at a relative humidity between 30-80%. Do not wet. Do not store near sterilizing agents.

Shelf Life

PCD-A chemical indicators have an expiration date of 5 years from the date of manufacture when stored under recommended storage conditions. Both, the expiration date and batch number of the product, are specified on the package label and in product unit.

Do not use PCD-A chemical indicators after their expiration date.

Do not reuse.

Endpoint stability reaction: PCD-A chemical indicator endpoint shall remain unchanged for a period of 6 months when stored under previously indicated conditions.

Disposal

Discard used chemical indicators along with paper waste, according to your country's healthcare and safety regulations.

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Terragene S.A.