Process challenge device (PCD)

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
128 °C \_\_\_\_\_ min.

135 °C \_\_\_\_\_ min.

SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

According to ISO 11140-1:2014 standard.

Lic. Adrián J. Rovetto  
 Director Técnico  
 Technical Director

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NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

Type 5 Integrator Indicator according to ISO 11140-1:2014

Indicador Integrador Tipo 5 conforme ISO 11140-1:2014

**Product description**

Process Challenge Device Integron® PCD26-2 Test Pack was developed to detect inadequate air removal and steam penetration in dynamic-air removal (pre-vacuum) steam sterilizers at 132/135 °C ≥ 4 minutes and gravity displacement steam sterilizers at 121 °C ≥ 30 minutes. The Test Pack consists of a porous barrier system equivalent to the challenge Test Pack specified in ANSI/AAMI ST79, holding in the middle a self-adhesive Record Card printed with Type 5 Integrator Indicator, where information about the sterilization cycle can be recorded. The indicator ink present on this self-adhesive card was developed to change to a reference color when a theoretical spore population reaches its time of death, indicating integration conditions were achieved. This condition is calibrated with the kill time of a 10<sup>6</sup> *Geobacillus stearothermophilus* ATCC 7953 spore population, calculated in BIER (Biological Indicator Evaluator Resistometer). It also has a Type 1 Process Indicator located on the outside of the package, which allows to quickly identify if it has been exposed to the sterilization process, and to distinguish between processed and unprocessed units.

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
**Unprocessed**

**Failures**

One or more parameters of sterilization were not achieved. The materials cannot be considered sterile and must be processed again. These results are usually obtained when steam quality is poor or adequate steam penetration or vacuum has not been generated.

**Processed and correct**