# PCD26-C

Rev. 11 / 12.2023

## Process Challenge Device

### Quality certification

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STEAM

LOT

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_


### IT26-C Chemical Integrator

LOT

\_\_\_\_\_

#### Stated Values (SV)

Temperature	Time
121 °C _____	_____ min.
128 °C _____	_____ min.
135 °C _____	_____ min.

SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

According to ISO 11140-1:2014 standard.

Terragene® and Bionova® are registered trademarks of Terragene S.A.  
ATCC® is a registered trademark of American Type Culture Collection.

Lic. Adhian J. Rovetto  
Director Técnico  
Technical Director

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#### Explanation of Symbols

**STEAM** Product designed for use in Steam sterilization cycles.  
**CE** CE mark.  
**EC REP** Authorized representative in the European Community.

**LOT** Batch number.  
 Manufacture Date.  
 Expiration Date.  
 Manufacturer.


**ATTENTION:** If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto** Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes.  
This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through **Bionova®** Software. Use **Trazanto** scanner to interpret the results (see instructions for use in the manual of the device).  
**ATTENTION:** do not write on the indicator before scanning.

#### Intended Use Table







Model	Autoclave/Steam Cycles
	Gravity Displacement
PCD26-C	121 °C
	30 minutes
	Dynamic Air Removal (Vacuum Assist)
	132 °C
	4 minutes

#### Result reference guide



##### Unprocessed

##### Accept

##### Reject

in a BIER (Biological Indicator Evaluator Resistometer). This pack also has a process indicator (Type 1 according to ISO 11140-1:2014 standard) printed on the outside of the package, which allows to quickly identify if it has been exposed to the sterilization process (color change from light blue to grey), and to distinguish between processed and unprocessed items.

#### Precautions

**WARNING:** Do not use PCD for monitoring Ethylene Oxide, Dry Heat, Formaldehyde or any sterilization process other than Steam.  
**WARNING:** Place one or more PCDs in sterilizing hard-to-reach areas to ensure all areas of the chamber are sterilized. Evaluate all load configurations to ensure ALL hard-to-reach areas have been identified, and place a PCD in each of those locations.  
**WARNING:** Do not use the sterilizer until the test result is correct.

#### Instructions for Use

- Place the pack inside a normally loaded steam autoclave, in those areas which are considered most inaccessible for the sterilizing agent (e.g., the center of the load and areas near the door).
- Run the sterilization cycle.
- After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the test pack. **NOTE:** The color of the box may vary from the original after undergoing the sterilization cycle. This does not represent a problem regarding the operation or quality of the product.
- Check that the process indicator printed on box has changed color from light blue to grey. Open the test pack and remove the integrating indicator. **PRECAUTION:** Wear safety glasses and gloves when removing the indicator from the sterilized test pack.
- Check IT26-C moving from chemical integrator for correct exposure. If the dark bar has reached the ACCEPT zone, this confirms that the inside of the pack has been exposed to sufficient sterilization conditions. For chemical integrator reference change, please refer to Result Reference Guide. Otherwise, check the sterilization process.
- Fill out the required information on the Record Card.
- Discard the pack immediately.

#### Monitoring frequency

Follow facility Policies and Procedures which should specify a chemical indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As the best practice and to provide optimal patient safety, Terragene® recommends that every sterilization load be monitored with an appropriate indicator.

#### Storage

Store in a dark place, at temperatures between 10-30 °C, 30-80 % Relative Humidity. Do not freeze. Do not store near sterilizing agents or other chemical products.

#### Shelf Life

PCD26-C has an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Do not use after expiration date.  
Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

#### Disposal

Discard the PCD after use along with paper waste, according to your country's healthcare and safety regulations.

## EN Process Challenge Device

For Steam sterilization

#### Indications for Use

##### United States

Terragene® Integron® PCD (PCD26-C) provides a defined challenge resistance against the claimed cycles shown below and demonstrated resistance equivalence to the AAMI/ANSI 16 towel pack. The device provides routine monitoring and sterilizer qualification testing steam sterilization processes. See Intended Use Table for reference.

##### Outside the United States

Integron® PCD26-C Steam Process Challenge Device has been designed for a rapid and easy monitoring of dynamic air removal Steam sterilization processes at 132/135 °C ≥ 4 minutes and gravity air displacement processes at 121 °C ≥ 30 minutes.

#### Device Description

Integron® PCD26-C Process Challenge Device (PCD) Test Pack has been designed to simulate a load to be sterilized and to pose a challenge to the sterilization process. It is used to evaluate the effective performance of the process by detecting inadequate air removal and steam penetration. Furthermore, it allows to perform the routine monitoring and periodic validation of the sterilizers (after repair, installation, relocation).  
Integron® PCD26-C consists of a disposable pre-assembled package as outlined in ANSI/AAMI ST79 which contains in the middle an IT26-C moving front chemical integrator (Type 5 according to ISO 11140-1:2014 standard) and a PCD26-C-RC self-adhesive Record Card where sterilization cycle information can be recorded. The extent of migration of the dark bar along the ACCEPT/REJECT zone of the integrating indicator indicates whether sterilization conditions were met or not. The ACCEPT result is reached when a theoretical spores population reaches its kill time, indicating integration condition has been reached. This condition is calibrated with the kill time of a 10<sup>6</sup> *Geobacillus stearothermophilus* ATCC® 7953 spores population, calculated

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