Infection Control Division | **Biological Indicators**

## BT110 Self-Contained Biological Indicator Rapid Readout Fluorescence System

### Usage

Monitoring Ethylene Oxide Sterilization Processes

### Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11138-1:2017 and ISO 11138-2:2017; IRAM 37102-1:1999 and IRAM 37102-2:1999.

### Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

### Classification

Class 1, according to risk.

### FDA 510(k)

K191021

### Characteristics

Polypropylene tube: 50.4 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm.

Polypropylene cap: 16.4 mm high x 10.7 mm external diameter. Wall thickness: 0.9 mm.

Cap filter: medical grade paper, 17.0 mm diameter.

Glass ampoule: 35.0 - 40.0 mm high. External diameter: 6.8 mm. Wall thickness: 0.2 - 0.3 mm.

Culture medium 0.5 - 0.7 ml, blue color.

Polypropylene microfibers on spore carrier, 17.0 mm diameter.

Spore carrier: filter paper, 16.0 mm diameter.

$\geq 10^6$  *Bacillus atrophaeus* ATCC 9372 spores per vial.

Final fluorescence reading is performed after 4 hour-incubation at 37 °C (sensitivity:  $\geq 97\%$ ).

An optional visual pH color change confirmation could be made after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at 37 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain blue after the incubation process.

7-day readout for visual confirmation by color change is optional and not intended to be routinely performed; it is an initial validation of the 4 hour-reading. Fluorescence results may be compared to the 7-day visual reading.

NOTE: If 7-day readout is performed, a humidified environment will be required to avoid medium dry out.

### Environmental conditions during manufacture

T= 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

Infection Control Division | **Biological Indicators**

## BT110 Self-Contained Biological Indicator Rapid Readout Fluorescence System

**Storage conditions**

T= 10-30 °C, RH 30-80 %, keep in the dark place in its original box.

**Transport conditions**

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent any risk for human health.

**Shelf life**

2 years.

**Packing**

50 units per box.

Packing information: product code and description, process for intended use, presentation, regulation, bacterial strain, storage conditions, manufacturer information and data on pack's label.

**Labelling**

On product: 17.0 mm x 33.0 mm polypropylene label. Chemical indicator line, printed with EO reactive ink (Color change to green). Graph showing final fluorescence reading time, product code, batch number, expiration date, process for intended use and name of organism printed in black.

On product's pack: product code and description, batch number, bacterial load, manufacture and expiration date, barcode and datamatrix code.

NOTE: manufacture date is calculated by subtracting 24 months to the expiration date.

**Possible target markets**

Healthcare and Industry.

**Other important information**

Incubation at 37 °C must be performed in Bionova® Auto-Reader Incubators. Read product's instructions for use thoroughly before use.

**Precautions**

Do not store the product near sterilizing agents.

Do not expose this product to Steam, Dry Heat, Radiation or any sterilization process other than EO.

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