Infection Control Division | Biological Indicators

BT20 Self-Contained Biological Indicator For steam sterilization processes

Usage

Monitoring vacuum-assisted and gravity air-displacement steam sterilization processes at 121-135 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11138-1:2017 and ISO 11138-3:2017; IRAM 37102-1:1999 and IRAM 37102-3:1999.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

Classification

Class 1, according to risk (ANMAT).

Characteristics

Polypropylene tube: 45.2 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm

Polypropylene cap: 13.0 mm high x 10.9 mm external diameter. Wall thickness: 0.7 mm

Cap filter: medical grade paper, 17.0 mm diameter.

Glass ampoule: 35.0 - 40.0 mm high. External diameter: 6.5 mm. Wall thickness: 0.2 - 0.3 mm.

Culture medium 0.5 - 0.7 ml, purple color.

Spore carrier: filter paper, 7.0 mm diameter.

10^5 or 10^6 *Geobacillus stearothermophilus* ATCC® 7953 spores per vial.

Upon completion of the sterilization cycle, incubate the indicator for a maximum of 24 hours at a temperature between 55-62 °C.

Positive growth results (failure in the sterilization process): culture medium turns to yellow after incubation between 55-62 °C for 24 hours (maximum).

D-value: not lower than 1.5 minutes at 121 °C. D-values at 132 °C and 135 °C are also informed.

Z-value: not lower than 6 °C.

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep in a dark place in its original box.

*ATCC® is a registered trademark of American Type Culture Collection.

Infection Control Division | **Biological Indicators**

BT20 Self-Contained Biological Indicator For steam sterilization processes

Transport conditions

It is recommended to store it away from sunlight and at a temperature between 10-30 °C.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

Acceptable excursions during transport are as follows:

T = 2-10 °C and T = 30-37 °C for no longer than seven days. No restrictions for relative humidity.

Shelf life

2 years.

Presentation

30 or 100 units per box.

10 units per bag.

Labelling

On product: 17.0 mm x 33.0 mm polypropylene label printed in black. 1.5 mm chemical indicator line, printed with steam reactive ink (color change to brown). Product code, batch number, expiration date, intended use, name of organism and blank spaces for data entry.

On product's box and label: product code and description, intended use, regulation, bacterial strain and load, storage conditions, manufacturer information, batch number, manufacture and expiration date, presentation, bar code and datamatrix code.

On product's bag: product code and description, intended use, regulation, bacterial strain and load, storage conditions, manufacturer information, batch number, manufacture and expiration date, presentation, bar code and datamatrix code.

Possible target markets

Healthcare and Industry.

Other important information

It is advisable to incubate at 60 °C in Bionova® IC10/20 incubator. Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, dry heat, vaporized hydrogen peroxide, radiation or any sterilization process other than steam.

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