Infection Control Division | **Biological Indicators**

## BT20S

Spore suspension for direct inoculation of products or preparation of customized biological indicators for monitoring Steam, Hydrogen Peroxide and Formaldehyde sterilization processes.

### Usage

For monitoring Steam, Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>) and Formaldehyde Sterilization Processes.

### Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11138-1:2017, ISO 11138-3:2017 and ISO 11138-5:2017; IRAM 37102-1:1999 and IRAM 37102-3:1999.

### Authorization

ANMAT (Argentinian National Administration of Drugs, Food and Medical Devices) PM 1614-1.

### Classification

Class 1, according to risk (ANMAT).

### Characteristics

Glass bottle: 18.0 mm diameter, 0.3 mm thick wall, 60.7 mm high.

Plastic cap: 16.0 mm external diameter, 9.4 mm high.

*Geobacillus stearothermophilus* ATCC 7953 spores.

D-Values<sub>STEAM</sub> are provided. Conditions: 121 °C, 132 °C and 135 °C.

Z-Values<sub>STEAM</sub> is provided.

D-Value<sub>FORM</sub> is provided. Conditions: 60 °C, 1 mol/l Formaldehyde.

D-Value<sub>H<sub>2</sub>O<sub>2</sub></sub> is provided. Conditions: 2 mg/l H<sub>2</sub>O<sub>2</sub>, 50 °C.

Upon completion of the sterilization cycle, the inoculated product should be tested in a suitable culture medium at 60 °C. We suggest the use of culture media Bionova® MC20-2 or MC1020-2 to final reading in 24 hours for Steam/H<sub>2</sub>O<sub>2</sub> sterilization and 48 hours for Formaldehyde sterilization. For conventional culture media like TSB, incubation for 7 days at 60 °C is recommended.

### Environmental conditions during manufacture

Refrigeration during inoculation. Sterility conditions are necessary only during the filling and sealing process, carried out under laminar flow conditions.

### Storage conditions

T = 2 - 8 °C (in refrigerator), RH 35-60 %. Keep out of direct light.

Infection Control Division | **Biological Indicators**

## BT20S

Spore suspension for direct inoculation of products or preparation of customized biological indicators for monitoring Steam, Hydrogen Peroxide and Formaldehyde sterilization processes.

### Transport conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent any risk for human health. Bottles can be stored at temperature below 30 °C for up to 15 days during transportation.

### Shelf life

2 years.

### Packing

1 bottle with 10 ml of spore suspension.

Packing information: product description, presentation, applicable regulation, manufacturer information and storage conditions.

### Labelling

On product's packing: product code and description, process for intended use, strain and population of bacteria, batch number, manufacture and expiration dates, barcode and datamatrix code.

On the bottle: Product code and description, presentation, bacterial strain and spore concentration, intended use, batch number, manufacture and expiration dates and manufacturer information.

### Possible target markets

Healthcare and Industry.

### Other important information

It is recommended to incubate at 60 °C in Bionova® IC10/20 incubator.

Read product's instructions for use thoroughly before use.

### Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam, H<sub>2</sub>O<sub>2</sub> and/or Formaldehyde.

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