# BT21/5

## Self-Contained spore ampoules.

*For Steam sterilization processes of liquids.*

### Usage

Monitoring industrial Steam sterilization processes of liquids between 121 °C - 135 °C.

### Applicable Regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11138-1:2017 and ISO 11138-3:2017; IRAM 37102-1:1999 and IRAM 37102-3:1999.

### Classification

Class 1, according to risk.

### Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

### Characteristics

Glass ampoule. 6.8 mm internal diameter. Wall thickness of 0.2 – 0.3 mm. Length: 40.0 mm.

Culture medium 0.5 ml, **purple** color.

10<sup>5</sup> *G. stearothermophilus* ATCC 7953 spores per ampoule.

Incubation conditions: 48 hours (maximum), between 55 - 62 °C.

Culture medium turns to **yellow** due to bacterial growth.

D-Value ≥ 1.5 minutes at 121 °C. Another D-value is declared at 132 °C and 135 °C.

Z-Value is informed.

### Environmental conditions during manufacture

T= 15-30 °C, RH= 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

### Storage conditions

T= 2-8 °C (in refrigerator), RH= 30-80 %, keep out of direct light. Storage at T ≤ 25 °C is possible for a period no longer than 15 days.

### Transport conditions

Storage conditions should be strictly followed.

Transport in closed reinforced boxes to avoid damage.

The transport of this product does not represent a risk for health.

Ampoules can be stored at temperature below 25 °C for up to 15 days during transportation.

### Shelf life

2 years.

### Packing

200 units + 20 negative controls per box (negative controls are marked with a colored dot on one end).

Packing information: product description, process for intended use, regulation, storage conditions, presentation, manufacturer information and data on box's label.

**Labelling**

On product's packing: product code and description, bacterial strain and load, batch number, manufacture and expiration date, bar code and datamatrix code.

**Possible target markets**

Healthcare and Industry.

**Other important information**

It is advisable to incubate at 60 °C in BIONOVA® IC10/20 Incubator.  
Read product's instructions for use thoroughly before use.

**Precautions**

Do not store the product near sterilizing agents.  
Do not use this product for monitoring EO, Dry Heat, Radiation or any sterilization process other than Steam.

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