BT21/6

Self-Contained spore ampoules.

For Steam sterilization processes of liquids.

Usage

Monitoring industrial Steam sterilization processes of liquids between 121 °C - 135 °C.

Applicable Regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11138-1:2017 and ISO 11138-3:2017; IRAM 37102-1:1999 and IRAM 37102-3:1999.

Classification

Class 1, according to risk.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

Characteristics

Glass ampoule. 6.8 mm internal diameter. Wall thickness of 0.2 – 0.3 mm. Length: 40.0 mm.

Culture medium 0.5 ml, **purple** color.

10⁶ *G. stearotherophilus* ATCC 7953 spores per ampoule.

Incubation conditions: 48 hours (maximum), between 55 - 62 °C.

Culture medium turns to **yellow** due to bacterial growth.

D-Value ≥ 1.5 minutes at 121 °C. Another D-value is declared at 132 °C and 135 °C.

Z-Value is informed.

Environmental conditions during manufacture

T= 15-30 °C, RH= 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

Storage conditions

T= 2-8 °C (in refrigerator), RH= 30-80 %, keep out of direct light. Storage at T ≤ 25 °C is possible for a period no longer than 15 days.

Transport conditions

Storage conditions should be strictly followed.

Transport in closed reinforced boxes to avoid damage.

The transport of this product does not represent a risk for health.

Ampoules can be stored at temperature below 25 °C for up to 15 days during transportation.

Shelf life

2 years.

Packing

200 units + 20 negative controls per box (negative controls are marked with a colored dot on one end).

Packing information: product description, process for intended use, regulation, storage conditions, presentation, manufacturer information and data on box's label.

Labelling

On product's packing: product code and description, bacterial strain and load, batch number, manufacture and expiration date, bar code and datamatrix code.

Possible target markets

Healthcare and Industry.

Other important information

It is advisable to incubate at 60 °C in BIONOVA® IC10/20 Incubator.
Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.
Do not use this product for monitoring EO, Dry Heat, Radiation or any sterilization process other than Steam.

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