Infection Control Division | **Biological Indicators**

BT225 Photon Biological Indicator

7-second readout fluorescence system

Usage

United States

Terragene® Bionova® Photon Biological Indicator (BT225) is a 7-second readout Self-Contained Biological Indicator (SCBI) inoculated with a minimum of 10^6 viable *Geobacillus stearothermophilus* bacterial spores and is intended for routine monitoring of the efficacy of steam sterilization processes. BT225 SCBI is not recommended and should not be used for qualification testing purposes. On each Terragene® Bionova® SCBI is a chemical process indicator that changes color from pink to brown when exposed to steam.

- Gravity-displacement Steam Sterilization Cycles
 - 132 °C, 25 minutes
 - 132 °C, 15 minutes
 - 132 °C, 10 minutes
 - 135 °C, 10 minutes

- Dynamic-air-removal Steam Sterilization Cycles
 - 132 °C, 4 minutes
 - 135 °C, 3 minutes

Terragene® Bionova® Photon Auto-reader Incubator (BPH) incubates at 60 °C and reads the Terragene® Bionova® Photon SCBIs at the times prescribed in the User Manual.

Outside the United States

Terragene® Bionova® Photon Biological Indicators (BT225) have been designed for quick and easy routine monitoring of dynamic-air-removal and gravity displacement steam sterilization cycles at 132-135 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN, ISO 13485:2016, ISO 11138-1:2017 and ISO 11138-3:2017; IRAM 37102-1:1999 and IRAM 37102-3:1999.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

Classification

Class 1, according to risk (ANMAT).

FDA De Novo

DEN220042.

Infection Control Division | Biological Indicators

BT225 Photon Biological Indicator

7-second readout fluorescence system

Characteristics

Polypropylene tube: 50.7 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm.

Polypropylene cap: 16.4 mm high x 10.7 mm external diameter. Wall thickness: 0.9 mm.

Cap filter: medical grade paper, 17.0 mm diameter.

Glass ampoule: 35.0 - 40.0 mm high. External diameter: 6.5 mm. Wall thickness: 0.2 - 0.3 mm.

Culture medium 0.4 - 0.6 ml, purple color.

Spore carrier: filter paper, 7.0 mm diameter.

≥ 10⁶ *Geobacillus stearothermophilus* ATCC®* 7953 spores per vial.

Incubate BT225 at (60±2) °C to get the final fluorescence result instantly (7-second readout) (sensitivity ≥ 97 %).

An optional visual pH color change confirmation could be made after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at 60 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain purple after the incubation process.

7-day readout is optional and not intended to be routinely performed; it is an initial validation of the instant readout. Fluorescence results may be compared to the 7-day visual reading. NOTE: if 7-day readout is performed, a humidified environment will be required to avoid medium to dry out.

D-Value: Not lower than 1.5 minutes at 121 °C. Not lower than 10 seconds at 132 °C. Not lower than 8 seconds at 135 °C.

Z-Value: Not lower than 10 °C.

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %. Sterility conditions are only necessary during the inoculation and manufacturing process, which is performed in laminar flow.

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep in a dark place in its original box.

Transport conditions

It is recommended to store it away from sunlight and at a temperature between 10-30 °C.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

Acceptable excursions during transport are as follows:

T = 2-10 °C and T = 30-37 °C for no longer than seven days. No restrictions for relative humidity.

*ATCC® is a registered trademark of American Type Culture Collection.

Infection Control Division | **Biological Indicators**

BT225 Photon Biological Indicator

7-second readout fluorescence system

Shelf life

18 months.

Packing

50 units per box.

Packing information: product code and description, process for intended use, presentation, regulation, bacterial strain, storage conditions, manufacturer information and data on pack's label.

Labelling

On product: 17.0 mm x 33.0 mm polypropylene label. 1.5 mm x 10.0 mm chemical indicator line, printed with steam reactive ink (color change to brown). Graph showing final fluorescence reading time, product code, batch number, expiration date, process for intended use and name of organism printed in black.

On product's packing: product code and description, batch number, bacterial load, manufacture and expiration date, barcode and datamatrix code.

Possible target markets

Healthcare and Industry.

Other important information

Incubation at 60 °C must be performed in Terragene® Bionova® Photon Auto-reader Incubator (BPH).

Read product's instructions for use thoroughly before use.

Terragene® Bionova® Photon Biological Indicator is intended to be marketed as Over-the-Counter (OTC) device.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to ethylene oxide, dry heat, radiation or any sterilization process other than steam.

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