Infection Control Division | **Biological Indicators**

BT225 Self-Contained Biological Indicator

Instant readout fluorescence system

Usage

Monitoring vacuum assisted and gravity air-displacement Steam sterilization cycles at 132-135 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN, ISO 13485:2016.

ISO 11138-1:2017 and ISO 11138-3:2017; IRAM 37102-1:1999 and IRAM 37102-3:1999.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

Classification

Class 1, according to risk (ANMAT).

Characteristics

Polypropylene tube: 50.7 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm.

Polypropylene cap: 16.4 mm high x 10.7 mm external diameter. Wall thickness: 0.9 mm.

Cap filter: medical grade paper, 17.0 mm diameter.

Glass ampoule: 35.0 - 40.0 mm high. External diameter: 6.5 mm. Wall thickness: 0.2 - 0.3 mm.

Culture medium 0.4 - 0.6 ml, purple color.

Spore carrier: filter paper, 7.0 mm diameter.

$\geq 10^6$ *Geobacillus stearothermophilus* ATCC®* 7953 spores per vial.

Incubate BT225 at (60±2) °C to get the final fluorescence result instantly (7 second readout) (sensitivity ≥ 97 %).

An optional visual pH color change confirmation could be made after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at 60 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain purple after the incubation process.

7-day readout is optional and not intended to be routinely performed; it is an initial validation of the instant readout.

Fluorescence results may be compared to the 7-day visual reading.

NOTE: if 7-day readout is performed, a humidified environment will be required to avoid medium to dry out.

D-Value: Not lower than 1.5 minutes at 121 °C. Another D-value is informed at 132 °C and at 135 °C.

Z-Value is informed.

*ATCC® is a registered trademark of American Type Culture Collection.

Infection Control Division | **Biological Indicators**

BT225 Self-Contained Biological Indicator

Instant readout fluorescence system

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %. Sterility conditions are only necessary during the inoculation and manufacturing process, which is performed in laminar flow.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place in its original box.

Transport conditions

It is recommended to respect the storage conditions during transport. Products should be transported in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent any risk for human health.

Shelf life

18 months.

Packing

50 units per box.

Packing information: product code and description, process for intended use, presentation, regulation, bacterial strain, storage conditions, manufacturer information and data on pack's label.

Labelling

On product: 17.0 mm x 33.0 mm polypropylene label. 1.5 mm x 10.0 mm chemical indicator line, printed with Steam reactive ink (Color change to brown). Graph showing final fluorescence reading time, product code, batch number, expiration date, process for intended use and name of organism printed in black.

On product's packing: product code and description, batch number, bacterial load, manufacture and expiration date, barcode and datamatrix code.

NOTE: manufacture date is calculated by subtracting 18 months to the expiration date.

Possible target markets

Healthcare and Industry.

Other important information

Incubation at 60 °C must be performed in Bionova® Auto-Reader incubators.

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to Ethylene Oxide, Dry Heat, Radiation or any sterilization process other than Steam.

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