BT30

Semi Self-contained Biological Indicator.

For Dry Heat sterilization processes.

Usage

Monitoring Dry Heat Sterilization Processes at 160 – 180 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11138-1:2017 and ISO 11138-4:2017; IRAM 37102-1:1999.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

Classification

Class 1, according to risk (ANMAT).

Characteristics

Tube A: nylon tube: 45.0 mm high x 8.4 mm external diameter. Wall thickness: 1.0 mm.

Tube B: Polypropylene random copolymer tube: 45.2 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm.

Cap A: nylon cap: 10.5 mm external diameter, 13.5 mm high. Wall thickness: 1.2 mm. This cap contains a 28.5 mm rod with a groove in which the filter paper is located.

Cap B: Polypropylene cap: 10.5 mm external diameter, 13.1 mm high. Wall thickness: 1.2 mm.

Filter paper: 5.0 x 25.0 mm, located in cap A rod.

Glass ampoule: 35.0 - 40.0 mm high. External diameter: 6.8 mm. Wall thickness: 0.2 – 0.3 mm. In tube B.

Culture medium: 0.5 – 0.7 ml, **blue** color.

10⁶ *Bacillus atrophaeus* ATCC 9372 spores per vial.

Positive growth results (failure in the sterilization process): Culture medium turns to **yellow**, incubation at 37 °C for 48 hours (maximum).

D-Value: Not lower than 2.0 minutes at 160 °C. Another D-Value is estimated at 170 °C for the calculation of Z-Value.

Z-Value: Not lower than 20 °C.

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep in a dark place.

Transport conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent any risk for human health.

Shelf life

24 months.

Packing

50 units per box - 25 tubes A + 25 tubes B.

Packing information: product code and description, process for intended use, presentation, classification according to regulation, manufacturer information and data on box's label. Note: manufacture date is calculated by subtracting 24 months to the expiration date.

Labelling

On product: 17.0 mm x 33.0 mm paper label. Printed in black. 1.5 mm chemical indicator line, printed with Dry Heat reactive ink (Color change to brown). Product code and batch number, process for intended use and bacterial strain.

On product's box: product code, batch number, bacterial load and manufacture and expiration date.

Possible target markets

Healthcare and Industry.

Other important information

It is advisable to incubate at 37 °C in BIONOVA® IC10/20 incubator. Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to Steam, Radiation or any sterilization process other than Dry Heat.

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