Infection Control Division | **Biological Indicators**

BT31

Spore ampoule for Dry Heat sterilization processes

Usage

Monitoring Dry Heat sterilization process. Compatible with extremely high temperatures generated during depyrogenation cycles and sterilization processes in ovens and dry heat tunnels.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.
ISO 11138-1:2017 and ISO 11138-4:2017; IRAM 37102-1:1999.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

Classification

Class 1, according to risk (ANMAT).

Characteristics

Glass ampoule. 10.2 mm internal diameter. Wall thickness of 0.5 mm. Length: 58.0 mm.

$\geq 10^6$ *Bacillus atrophaeus* ATCC 9372 spores per ampoule.

Upon completion of the sterilization cycle, the content of the ampoule should be incubated in a suitable culture medium at 37 °C. We suggest the use of Bionova® MCI030-2 culture medium to final reading in 72 hours.

D-Value_{DRY} is provided: Not lower than 2.0 min. Conditions: 160 °C.

Z-Value_{DRY} is provided: Not lower than 20 °C

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

Storage conditions

T = 10 – 30 °C, RH 30-80 %, keep in a dark place.

Transport conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent any risk for human health.

Shelf life

2 years.

Packing

50 units per box.

Packing information: product code and description, process for intended use, regulation, storage conditions, presentation, manufacturer information and data on box's label.

Infection Control Division | **Biological Indicators**

BT31

Spore ampoule for Dry Heat sterilization processes

**Labelling**

On the ampoule: product code, process for intended use, name of organism, batch number, expiration date and trademark.

On product's packing: product code, bacterial load, batch number, expiration date, barcode and datamatrix code.

Possible target markets

Healthcare and Industry.

Other important information

It is advisable to incubate the content of BT31 in MC1030-2 culture medium at 37 °C using Bionova® IC10/20 incubators.

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to Steam, Radiation or any sterilization process other than the recommended one.

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