# BT70S/5

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# BT70S/5

## Spore suspension

For direct inoculation of products or preparation of customized biological indicators for monitoring Gamma Radiation Sterilization Processes.

## Usage

For monitoring Gamma Radiation Sterilization Processes. For industrial use only.

## Applicable Regulation

ISO 11138-1:2006; IRAM 37102-1:1999.

## Classification

Class 1, according to risk.

## Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

## Characteristics

Glass bottle: 18 mm diameter, 0.3 mm thick wall, 60.7 mm high.

Plastic cap: 16 mm external diameter, 9.4 mm high.

Population:  $10^4$  *Bacillus pumilus* ATCC 27142 spores per 0.1 ml ( $10^5$  spores per ml), prepared in 40% ethanol.

D-Value is provided. Calculated using Limited Holcomb-Spearman-Karber method.

Upon completion of the sterilization cycle, the inoculated product should be tested in a suitable culture medium at 37 °C. We suggest the use of culture media Bionova® MC70 to final reading in 48 hours. For conventional culture media like TSB, incubation for 7 days at 37 °C is recommended.

## Environmental conditions during manufacture

Refrigeration during inoculation. Sterility conditions are necessary only during the filling and sealing process, carried out under laminar flow conditions.

## Storage conditions

For short-term use, store at 2–8 °C in refrigerator. Freeze to store for long periods (-10 to -25 °C). RH: 35-60%. Keep out of direct light.

## Transportation conditions

Storage conditions should be strictly followed.

Transport in closed reinforced boxes to avoid damage. This product does not represent a risk for health. Bottles can be stored at temperature below 30 °C for up to 15 days during transportation.

## Shelf-life

24 months.

**Packing**

1 bottle with 10 ml of spore suspension within a plastic bag.

**Labelling**

On the bottle: Product code and description, presentation and spore concentration, directions for use, batch number, manufacture and expiration date and manufacturer information.

**Possible target markets**

Healthcare and Industry

**Other important information (optional)**

It is recommended to incubate at 37 °C in BIONOVA® IC10/20 incubator.

Read product's directions for use thoroughly before use.

**Note:** When necessary, limits of each product's parameters should be specified.

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to Steam, H<sub>2</sub>O<sub>2</sub>, EO, Dry Heat or any sterilization process other than Radiation.

The positive biological indicator can be sterilized in gravity air displacement steam sterilizers at 121°C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.

Note: Spores can be incubated in culture media BIONOVA® MC70 culture medium or other culture medium made up of soybean casein and yeast extract.