Infection Control Division | **Biological Indicators**

BT98 Self-Contained Biological Indicator Hyper Rapid Readout Fluorescence System

Usage

United States

Bionova® Hyper Biological Indicator (BT98) is a Self-Contained Biological Indicator (SCBI) inoculated with a minimum of 10⁶ viable *Geobacillus stearothermophilus* bacterial spores and is intended for monitoring the efficacy of vaporized hydrogen peroxide sterilization processes in the following systems:

Vaporized hydrogen peroxide	
Cycles	
Sterrad™ 100S Sterilization System	
Sterrad NX™ Sterilization System (Standard and Advanced Cycles)	
Sterrad NX™ with ALLClear™ Technology Sterilization System (Standard and Advanced Cycles)	
Sterrad™ 100NX Sterilization System (Standard, Flex, Express and Duo Cycles)	
Sterrad™ 100NX with ALLClear™ Technology Sterilization System (Standard, Flex, Express and Duo Cycles)	
V-Pro® S2 Low Temperature Sterilization System (Fast, Non Lumen, Lumen and Flexible Cycles)	
V-Pro® maX 2 Low Temperature Sterilization System (Fast Non Lumen, Non Lumen, Lumen and Flexible Cycles)	
Fluorescence read time	pH color change
5 minutes	5 days

Bionova® BT98 has hyper rapid readout at 5 minutes at 60 °C.

Outside the United States

Bionova® BT98 fluorescence hyper rapid readout biological indicators have been designed for quick and easy monitoring of vaporized hydrogen peroxide sterilization cycles.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/EN ISO 13485:2016. ISO 11138-1:2017 and ISO 11140-1:2014.

Authorization

ANMAT (Argentinian National Administration of Drugs, Food and Medical Devices) PM 1614-1.

Classification

Class 1, according to risk (ANMAT).

FDA 510(k)

K221641.

Infection Control Division | Biological Indicators

BT98 Self-Contained Biological Indicator Hyper Rapid Readout Fluorescence System

Characteristics

Polypropylene tube: 50.4 mm high x 8.5 mm external diameter.
Wall thickness: 0.5 mm.

Polypropylene cap: 16.4 mm high x 10.7 mm external diameter. Wall thickness: 0.9 mm.

Cap filter: polyethylene fibers. 17.0 mm diameter, 64.5 g/m².

Glass ampoule: 35.0-40.0 mm high. External diameter: 6.8 mm.
Wall thickness: 0.2-0.3 mm.

Culture medium 0.5-0.7 ml, purple color.

Polypropylene microfibers on spore carrier: 17 mm of diameter, 60 g/m².

Spore carrier: polyethylene fibers. 17 mm of diameter, 64.5 g/m².

≥ 10⁶ *G. stearothermophilus* ATCC®* 7953 spores per vial.

Each BT98 has a process indicator on the label that changes from purple to green when exposed to hydrogen peroxide. The SCBI also has a Data Matrix on its label for unambiguous identification (compatible with Bionova® Q).

Final fluorescence readout is performed after 5-minute incubation at 60 °C (sensitivity ≥ 97 %).

An optional visual confirmation could be made after 5 days of incubation. Color change of the culture medium is developed by pH variation when microorganisms grow. If sterilization process has not been successful, culture medium will change from purple to yellow during incubation at 60 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain purple after the incubation process.

NOTE: If 5-day readout is performed, a humidified environment will be required to avoid medium to dry out.

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %. Sterility conditions are only necessary during the inoculation and manufacturing process, which is performed in laminar flow.

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep away from sunlight in its original box.

Transport conditions

It is recommended to store it away from sunlight and at a temperature between 10-30 °C.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

*ATCC® is a registered trademark of American Type Culture Collection.

Infection Control Division | Biological Indicators

BT98 Self-Contained Biological Indicator Hyper Rapid Readout Fluorescence System

Shelf life

2 years.

Packing

50 units per box.

Packing information: Product code and description, process for intended use, presentation, regulatory and standards information, bacterial strain, storage conditions and manufacturer information.

Labelling

On product's packing: Product code and description, batch number, bacterial load, manufacture and expiration date, barcode and Data Matrix.

Possible target markets

Healthcare and Industry.

Other important information

It is advisable to incubate at 60 °C in Bionova® Hyper (BHY) Auto-reader incubators.

Read product's Instructions for Use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not use this product for monitoring steam, ethylene oxide, dry heat, formaldehyde or other sterilization processes different from the specified.

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Peróxido de hidrógeno vaporizado	
Ciclos	
Sterrad™ 100S Sterilization System	
Sterrad NX™ Sterilization System (Standard and Advanced Cycles)	
Sterrad NX™ with ALLClear™ Technology Sterilization System (Standard and Advanced Cycles)	
Sterrad™ 100NX Sterilization System (Standard, Flex, Express and Duo Cycles)	
Sterrad™ 100NX with ALLClear™ Technology Sterilization System (Standard, Flex, Express and Duo Cycles)	
V-Pro® S2 Low Temperature Sterilization System (Fast, Non Lumen, Lumen and Flexible Cycles)	
V-Pro® maX 2 Low Temperature Sterilization System (Fast Non Lumen, Non Lumen, Lumen and Flexible Cycles)	
Tiempo de lectura por fluorescencia	Cambio de color por pH
5 minutos	5 días

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