# BTC90-2/6

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# BTC90-2/6

## Spore coupon

*For Plasma or Vaporized Hydrogen Peroxide*

### Usage

Monitoring plasma or vaporized Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>) sterilization processes. To be used in Helix-PCD systems.

### Applicable Regulation

ISO 11138-1:2006 and IRAM 37102-1:1999.

### Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

### Characteristics

- Stainless steel coupon of 38.50 mm x 6.0 mm x 0.8 mm
  - 10<sup>6</sup> *Geobacillusstearothermophilus* ATCC 7953 spores per coupon.
- Upon completion of the sterilization cycle, the spore coupon should be incubated in a suitable culture medium at 60 °C. We suggest the use of BIONOVA® MC20 or MC1020 culture media to final readout at 24 hours. For conventional culture media like TSB, it is recommended a 7-day incubation at 60 °C.
- D Value: at 50 °C, 2mg/l H<sub>2</sub>O<sub>2</sub>.

### Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

### Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place.

### Transportation conditions

Storage conditions should be strictly followed. Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

### Shelf-life

24 months.

**Packing**

100 units per box.

Packing information: Code, product description, process in which is used, presentation, classification according to regulation, batch number, production date, expiration date and manufacturer information.

Weight per box: 152.5g

**Labelling**

On the product's box: batch number, production date and expiration date.

**Possible target markets**

Healthcare and Industry.

**Other important information (optional)**

It is advisable to incubate at 60 °C in BIONOVA® IC10/20 incubator. Read product's directions for use thoroughly before use.

NOTE: When necessary, limits of each product's parameters should be specified.

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry heat, Radiation or any sterilization process other than Hydrogen Peroxide.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilizer, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.

Note: The product can be incubated in BIONOVA® MC20 or MC1020 culture media or in a culture medium made of soybean casein and yeast extract.