# KMB222

## Kit lectora incubadora Bionova® MiniBio + Indicadores Biológicos Auto-contenidos Bionova® BT222.

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# KMB222

## Bionova® BT222 Self-Contained Biological Indicators + Bionova® MiniBio auto-reader incubator Kit.

*For incubation and reading of Super Rapid Readout Biological Indicators Bionova® BT222*

### Usage

**BT222 Super Rapid Biological Indicators with Super Rapid Fluorescence Reading System:** control of steam sterilization processes by dynamic air removal (pre-vacuum) and steam sterilization by gravity at 121 °C - 135 °C.

**Bionova® Minibio Incubator reader:** reading and incubation of biological indicators Bionova® BT222, 1 hour at 60°C.

Additionally, it also allows the incubation and reading of the following indicators:

Bionova® Biological Indicators with Rapid Reading System by Fluorescence: BT220 (Steam): 3 hs at 60 °C, BT221 (Steam): 1 h. at 60 °C, BT110 (OE): 4 h. at 37 °C, BT102 (Formaldehyde): 2 h at 60 °C, BT95 (VH2O2): 2 h. at 60 °C.

Bionova® Biological Indicators with Super Rapid Fluorescence Reading System: BT223 (Steam): 30 min. at 60 °C, BT96 (VH2O2): 30 min. at 60 °C.

Bionova® Biological Indicators with Ultra Rapid Reading System by Fluorescence: BT224 (Steam): 20 min. at 60°C

### Applicable regulation

MiniBio: Low Voltage Directive 2014/35/EU, EMC Directive 2014/30/CE and RoHS Directive 2011/65/EU.

BT222: ISO 11138-1: 2006 and ISO 11138-3: 2006; IRAM 37102-1: 1999 and IRAM 37102-3: 1999

### Authorization

ISO 13485:2003/NS-EN ISO 13485:2012. Bionova® BT222: ANMAT PM 1614-1.

### Classification

Bionova® BT222: Class 1, according to risk.

### Characteristics

#### MiniBio

3 positions (38 mm depth, 9 mm diameter) to incubate Rapid, Super Rapid and Ultra Rapid Self-contained Biological Indicators.

1 ampoule crusher.

Orifice for external thermometer.

Built-in top cover.

Dual Temperature System: allows to select between two incubation temperatures (37 °C and 60 °C)

USB connection for recording results on PC using a Reading and Traceability software.

It allows to incubate 3 indicators simultaneously, in different programs, at the same temperature.

Thermal Printer of Results.

Audible alarm for indication of events

Detection and automatic cancellation of the fluorescence reading of the biological indicator.

Fluorescence reading wavelength: ~ 460 nm

Firmware update by the user.

Dimensions: 180 mm high, 154 mm wide and 105 mm deep.

Voltage Range: 100 - 240 V AC

Power: 28 W

Frequency: 50 - 60 Hz

**BT222**

Polypropylene tube: 50.4 mm high x 8.5 mm external diameter. Wall of 0.5 mm thick.

Polypropylene microfiber cap on the base of the tube, 17.0 mm in diameter, 80 g / m<sup>2</sup>.

Polypropylene plastic lid: 10.5 mm external diameter, 16.4 mm high. Wall of 0.9 mm thick

Glass ampoule: 35.0 to 40.0 mm in height. External diameter: 6.8 mm. Culture medium 0.5 - 0.7 ml, purple color.

Filter paper: medical grade paper 60 g / m<sup>2</sup>, diameter 17.0 mm.

Spore carrier: Paper filter: 5.0 x 25.0 mm, 160 g/m<sup>2</sup>.

≥ 10<sup>6</sup> spores of *Geobacillus stearothermophilus* ATCC 7953 per vial.

The final fluorescence reading is carried out after 1 hour of incubation at 60 °C (sensitivity: 97%).

Optionally, a visual confirmation can be made by changing the color by changing the pH after a 48-hour incubation. If the sterilization process has not been successful, the culture medium will change to a greenish color first, and then to yellow during the incubation at 60°C, indicating the presence of live spores. If the sterilization was successful, the culture medium will remain purple after the incubation process.

The reading at 7 days for confirmation is optional and does not need to be done routinely; it is an initial validation of the reading at 1 hour. The fluorescence results can be compared with the 7-day reading.

NOTE: if the reading is made after 7 days, a humidified environment will be required to prevent the medium from drying out.

Value D: not less than 1.5 minutes at 121 °C. Another D value is declared at 135°C

Z value: not lower than 6 °C.

**Environmental conditions during manufacture**

15-30 °C (room temperature), RH 30-80 %.

BT222: Conditions of sterility only during the inoculation process that is carried out under laminar flow.

**Storage conditions**

T = 10-30 °C, RH 30-80 %, store in dark place.

**Transportation conditions**

Storage conditions should be strictly followed. Transport in closed and strengthened boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

**Shelf-life**

MiniBio: 10 years. BT222: 2 years.

**Packing**

1 MiniBio reader-incubator + 12 aluminum bags with hermetic closure, containing 10 units of BT222 each.

Data on the packaging: code and description of the product, process for which it is used, presentation, manufacturer information and information on the label.

## Labelling

In the package of the kit: code, batch and description of the kit, code and batch of MiniBio and code and batch of BT222.

In the MiniBio box: same label on the product and on the packaging: product code and description, batch, serial number and manufacturer's data.

In the BT222 bag: product code, batch, bacterial population, date of manufacture and expiration.

In each vial of BT222: polypropylene label of 17.0 mm x 33.0 mm. Printed in black and purple colors. Chemical indicator line of 1.8 mm x 10.0 mm printed with Steam reactive ink (turns brown). Graph showing the final reading time by fluorescence. Product code, batch, process for which it is used and bacterial strain.

## Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

## Other important information (optional)

Read product's directions for use thoroughly before use

The Bionova® MiniBio has a button on the incubator board that allows you to select an incubation temperature of 37 °C or 60 °C. In addition, it has a button for each of the reading positions that allows you to select a different time program for each position. The Bionova® MiniBio incubator has an automatic reading system and a printer to record the data.

## Description of identified and/or predictable risks which could represent a risk for the product

### MiniBio

Do not pour any liquid inside.

Do not immerse the incubator in any liquid.

Use only indoors.

Disconnect the power cable before cleaning.

Do not use abrasive or corrosive cleaners or disinfectants.

Make sure the incubator is connected to a power cord with the proper power.

### BT222

Do not store the product near sterilizing agents.

Do not use the Biological Indicator to control sterilization processes by EO, Dry Heat, Radiation or any other process different from Steam sterilization.

Once the product is used it is recommended to sterilize in a steam sterilizer by gravity for a minimum of 20 minutes at 121°C or for a minimum of 15 minutes at 132°C, or in a vacuum steam sterilizer for a minimum of 10 minutes at 134 °C.