MC1020-2

Culture medium with pH indicator.

For growth of *Geobacillus stearothermophilus* ATCC 7953, *Bacillus subtilis* ATCC 35021, *Bacillus atrophaeus* ATCC 9372 and *Bacillus pumilus* ATCC 27142 bacterial spores.

Usage

To test growth of bacterial spores from the following strains: *Geobacillus stearothermophilus* ATCC 7953 used for monitoring Steam, Plasma or Vaporized Hydrogen Peroxide and Formaldehyde sterilization processes.

Bacillus subtilis ATCC 35021 used for monitoring Steam at low temperature sterilization processes.

Bacillus atrophaeus ATCC 9372 used for monitoring Dry Heat and EO sterilization processes.

Bacillus pumilus ATCC 27142 used for monitoring Gamma Radiation sterilization processes.

Applicable Regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11138-1:2017, ISO 11138-2:2017, ISO 11138-3:2017, ISO 11138-4:2017 and ISO 11138-5:2017; IRAM 37102-1:1999, IRAM 37102-2:1999 and IRAM 37102-3: 1999.

Classification

Class 1, according to risk.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

Characteristics

Sterile growth medium is specially formulated for use with Bionova® products (also suitable for other brand products).

Growth medium color: red. Volume: 2.0 ml

Plastic tube: 76.0 mm high, 12.5 mm diameter

Incubation conditions:

24 hours (maximum), between 55-62 °C for *Geobacillus stearothermophilus* ATCC 7953 spores subjected to Steam/VH₂O₂ processes and 48 hours (maximum) for Formaldehyde processes.

48 hours (maximum), between 35-39 °C for *Bacillus atrophaeus* ATCC 9372 spores subjected to Dry Heat and EO processes, *Bacillus subtilis* ATCC 35021 spores subjected to Steam at low temperature processes and *Bacillus pumilus* ATCC 27142 spores subjected to Gamma Radiation processes.

Medium turns to yellow when spores grow.

Environmental conditions during manufacture

T= 15-30 °C, RH= 30-80 %. Sterility conditions are necessary during plastic tube filling with sterile medium, performed in laminar flow.

Storage conditions

T= 4-20 °C, RH= 30-80 %, keep out of direct light.

Transport conditions

Storage conditions should be strictly followed.

Transport in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent a risk for health.

Shelf life

2 years.

Packing

50 units per box.

Packing information: product description, storage conditions, presentation, regulation, manufacturer information and data on box's label.

Labelling

On product: Label 17.0 mm x 42.0 mm printed in black. Code and description of the product, intended use, batch number and expiration date.

On product's packing: code and description of the product, intended use, color change, presentation, batch number, manufacture date, expiration date, bar code and data matrix code.

Possible target markets

Healthcare and Industry.

Other important information

It is advisable to incubate at 37/60 °C in BIONOVA® IC10/20 incubator.

Precautions

Do not store the product near sterilizing agents.

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