KPCD220-1

Steam Test Pack PCD Kit

Process Challenge Device for Steam sterilization processes

Usage

For routine load release, especially for loads containing implants. For routine sterilizer monitoring. For periodic validation of the sterilizer (after installation, repair, relocation).

Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity displacement steam autoclaves.

Condition: vacuum-assisted steam sterilization processes at 132/135 °C ≥ 4 minutes and gravity-air displacement steam sterilization processes at 121 °C ≥ 30 minutes.

Applicable Regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11140-1:2014, ISO 11138-1:2017, ISO 11138-3:2017, ANSI/AAMI ST79.

Classification

Class 1, according to risk.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Characteristics

Porous cards system holding a Bionova® BT220 Self-contained Biological Indicator and a Chemical Integrator (Type 5 according to ISO 11140-1: 2014) printed on PCDBI-1-RC self-adhesive Record Card, where data on sterilization cycle may be written. Two possibilities for documentation record: the whole Record Card can be adhered, or alternatively, only the area containing the indicator may be recorded. Cardboard box (10.5 cm x 6.5 cm x 1.8 cm approximately) with a Process Indicator (Type 1 according to ISO 11140-1: 2014).

Process Indicator:

Initial color: **blue**

Final color: **dark grey or black**

PCDBI-1-RC Chemical Integrator (within the challenge device):

Initial color: **yellow**

Final color: **black**

100 % Toxic Heavy Metals free

Integrating condition is calibrated against the kill time of a 10⁶ G. *stearothermophilus* ATCC 7953 spores' population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

BT220 Rapid Self-contained Biological Indicator (within the challenge device): ≥ 10⁶ G. *stearothermophilus* ATCC 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

Final results: fluorescence readout after a 3-hour incubation at 60 °C (sensitivity \geq 97%). A visual confirmation of the result is optional (color change of the culture medium caused by pH change) and can be made after incubation at 60 °C for 48 hours and/or 7 days.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

Environmental conditions during manufacture

T= 15-30 °C, RH= 30-80 %.

Storage conditions

T= 10-30 °C, RH= 30-80 %, keep out of direct light.

Transportation conditions

Storage conditions should be strictly followed. Carry in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

2 years.

Packing

25 PCD220-1 + 25 BT220 Biological Indicators per box.

Packing information: product description, storage conditions and manufacturer information.

Labelling

On product's packing: product code and description, presentation, process for intended used, regulation, batch number, batch number of the Self-Contained Biological Indicator, manufacture and expiration date, bar code and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information

Read product's instructions for use thoroughly before use.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

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