# PCD222-1

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# PCD222-1

## Process Challenge Device

Test Pack for Steam sterilization processes

### Usage

For routine load release, especially for loads containing implants. For routine sterilizer monitoring. For periodic validation of the sterilizer (after installation, repair, relocation)

Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity displacement steam autoclaves.

**Condition:** 132/135 °C ≥ 4 minutes and 121 °C ≥ 30 minutes, steam sterilization processes.

### Applicable Regulation

ISO 11140-1:2014, ISO 11138-1:2006, ISO 11138-3:2006, ANSI/AAMI ST79

### Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-4.

### Characteristics

Porous cards system holding a Self-contained Biological Indicator and a Type 5 Integrator Indicator printed on PCDBI-1-RC self-adhesive Record Card, where data on sterilization cycle may be written. Two possibilities for documentation record: the whole Record Card can be adhered, or alternatively, only the area containing the indicator may be recorded.

Cardboard box (10.5 cm x 6.6 cm x 1.8 cm) with a Type 1 Process Indicator printed.

### Process Indicator:

Initial color: **blue**

Final color: **dark grey or black**

### PCDBI-1-RC Integrator Indicator (within the challenge device):

Initial color: **yellow**

Final color: **black**

Integration condition is calibrated against the kill time of a  $10^6$  *G. stearothermophilus* ATCC 7953 spores' population, calculated in BIER (Biological Indicator Evaluator). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

**Self-contained Super Rapid Biological Indicator** (within the challenge device):  $\geq 10^6$  *G. stearothermophilus* ATCC 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

**Final results:** by fluorescence after incubation for 1 hour at 60 °C (sensitivity 97%). A 7-day readout is optional to confirm 1-hour results by color change.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

**Environmental conditions during manufacture**

15-30 °C (room temperature), RH 30-80 %.

**Storage conditions**

Best stored in original box under the next conditions: T = 10-30 °C, RH 30-80 %.

**Transportation conditions**

Storage conditions should be strictly followed. Carry in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

**Shelf-life**

2 years.

**Packing**

Packing information: product code and description, operating instructions, classification according to regulation, process for intended use, batch number, expiration date and manufacturer information.

Weight per unit: 217,27 g.

**Labelling**

Information only on product's box.

**Possible target markets**

Healthcare, Food, Pharmaceutical and Medical Industries.

**Other important information (optional)**

Read product's directions for use thoroughly before use.

**Note:** When necessary, limits for each product's parameters should be specified.

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.