# KPCD222-2

## Steam Test Pack PCD Kit

*Process Challenge Device for Steam sterilization processes*

### Usage

For routine load release, especially for loads containing implants. For routine sterilizer monitoring. For periodic validation of the sterilizer (after installation, repair, relocation).

Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity displacement steam autoclaves between 121-135 °C.

### Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11140-1:2014, ISO 11138-1:2017, ISO 11138-3:2017, ANSI/AAMI ST79.

### Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

### Classification

Class 1, according to risk (ANAMAT).

### Characteristics

The kit consist of:

- PCD222-2: Porous cards system holding:

A Self-contained Biological Indicator Bionova® BT222.

A Type 5 chemical indicator.

A self-adhesive Record Card, PCDBI-2-RC, where data on sterilization cycle may be written.

Two possibilities for documentation record: the whole Record Card can be adhered, or alternatively, only the area containing the indicator may be recorded.

- Self-contained Biological Indicator Bionova® BT222 (to be used as a positive control)

#### Process Indicator:

Initial color: **blue**.


Final color: **dark grey or black**.

#### PCDBI-2-RC Integrator Indicator (within the challenge device):

Initial color: **purple**.

Final color: **green**.

100 % Toxic Heavy Metals free.

Integration condition is calibrated against the kill time of a  $10^6$  *G. stearothermophilus* ATCC 7953 spore population, calculated in BIER (Biological Indicator Evaluator). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

**Self-contained Super Rapid Biological Indicator BT222** (within the challenge device):  $\geq 10^6$  *G. stearothermophilus* ATCC 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

**Final results:** by fluorescence after incubation for 1 hour at 60 °C (sensitivity 97%). A 7-day readout is optional to confirm 1-hour results by color change.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

### Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

### Storage conditions

T = 10-30 °C, RH = 30-80 %, keep out of direct light.

### Transport conditions

Storage conditions should be strictly followed. Carry in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent a risk for health.

### Shelf life

2 years.

### Packing

Presentations: 25 PCD222-2 + 25 BT222 Biological Indicators per box.

25 PCD222-2 + 5 BT222 Biological Indicators per box.

16 PCD222-2 + 4 BT222 Biological Indicators per box.

Packing information: product description, storage conditions, manufacturer information and data on box's label.

### Labelling

On product's packing: product code and description, presentation, regulation, batch number, batch number of the Self-contained Biological Indicator, manufacture and expiration date, bar code and datamatrix code.

### Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

### Other important information

Read product's instructions for use thoroughly before use.

### **Precautions**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.

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