Infection Control Division | **Chemical Indicators**

PCD224-2 Process Challenge Device.

Test pack for Steam sterilization processes.

Usage

For routine load release, especially for loads containing implants. For routine sterilizer monitoring. For periodic validation of the sterilizer (after installation, repair, relocation).

Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity displacement steam sterilization processes between 132-135 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11140-1:2014, ISO 11138-1:2017, ISO 11138-3:2017, ANSI/AAMI ST79.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Classification

Class 1, according to the risk (ANMAT).

FDA 510(k)

K191021

Characteristics

Porous cards system holding a Bionova® BT224 Self-contained Biological Indicator and a Chemical Integrator (Type 5 Integrating Indicator according to ISO 11140-1: 2014) printed on PCDBI-2-RC self-adhesive Record Card, where data on sterilization cycle may be written. Two possibilities for documentation record: the whole Record Card can be adhered, or alternatively, only the area containing the indicator may be recorded.

Cardboard box (10.5 cm x 6.5 cm x 1.8 cm approximately) with a Process Indicator printed (Type 1 according to ISO 11140-1: 2014).

Process Indicator

Initial color: **light blue**

Final color: **dark grey or black**

PCDBI-2-RC Chemical Integrator (within the challenge device)

Initial color: **purple**

Final color: **green**

100 % Toxic Heavy Metals free

Integration condition is calibrated against the kill time of 10^6 G. *stearothermophilus* ATCC 7953 spores' population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

Infection Control Division | **Chemical Indicators**

PCD224-2 Process Challenge Device.

Test pack for Steam sterilization processes.

BT224 Ultra Rapid Self-contained Biological Indicator (within the challenge device): $\geq 10^6$ *G. stearothermophilus* ATCC 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

Final results: fluorescence readout after 20-minutes incubation at 60 °C (sensitivity $\geq 97\%$). A visual confirmation of the result is optional (color change of the medium caused by pH change) and can be made after incubation at 60 °C for 48 hours and/or for 7 days.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep in a dark place.

Transport conditions

Storage conditions should be strictly followed. Products should be transported in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent any risk for human health.

Shelf life

2 years.

Packing

Packing information: product code and description, regulation, process for intended use, conditions and instructions for use, storage conditions, presentation, manufacturer data and label information.

Labelling

On the package: product code, batch number, batch number of the Self-Contained Biological Indicator, manufacture date, expiration date, regulation, bar code and datamatrix.

On the BT224 Self-Contained Biological Indicator: label printed with Steam reactive ink (turns to brown), product code, batch number, expiration date, process for intended use and name of the organism. On the Chemical Integrator indicator PCDBI-2-RC: product code, process for intended use, classification according to the regulation, color change, batch number and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Infection Control Division | **Chemical Indicators**

PCD224-2 Process Challenge Device.

Test pack for Steam sterilization processes.

Other important information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

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